

# Zimmer<sup>®</sup> Collagen Plug, Tape and Patch

## Protect. Heal. Repair.



### ABSORBABLE WOUND DRESSINGS

**1** Protects Wound Bed  
Adheres and provides coverage to oral wounds and sores

**2** Aids Wound Healing  
Porous, absorbable matrix supports delicate new tissue

**3** Highly Porous,  
Highly Absorbent  
Zimmer Collagen Plug holds up to 30x own weight in fluid  
Resorbs in fewer than 30 days  
Greater than 90% open pores

**4** Versatile For  
Everyday Use  
Three convenient shapes for common oral wounds and procedures

# ABSORBABLE COLLAGEN MATRICES

## IN THREE CONVENIENT SHAPES

Zimmer  
Collagen Plug



10 mm x 20 mm  
Cat. No. 0102Z

Zimmer  
Collagen Tape



2.5 cm x 7.5 cm, 1 mm thick  
Cat. No. 0100Z

Zimmer  
Collagen Patch



2 cm x 4 cm, 3 mm thick  
Cat. No. 0101Z

## VERSATILE WOUND PROTECTION

- Designed to resorb within a short timeframe
- Indicated for management of oral wounds and sores
  - Denture sores
  - Oral ulcers (non-infected or viral)
  - Periodontal surgical wounds
  - Suture sites
  - Burns
  - Extraction sites
  - Surgical wounds
  - Traumatic wounds



Atraumatic tooth extraction.



Modified Zimmer Collagen Plug delivered to protect grafted site.



Sutured site.



Ovate pontic contours soft tissue during healing.

Contact us at 1-800-342-5454 or visit [zimmerbiometdental.com](http://zimmerbiometdental.com)

Zimmer Biomet Dental  
Global Headquarters  
4555 Riverside Drive  
Palm Beach Gardens, FL 33410  
Phone: +1-561-776-6700  
Fax: +1-561-776-1272

All references to Zimmer Biomet Dental contained herein refer to the dental subsidiaries of Zimmer Biomet Holdings, Inc. Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured, distributed and marketed by Zimmer Biomet Dental (and, in the case of distribution and marketing, its authorized marketing partners). Zimmer Collagen Plug, Tape and Patch are manufactured by Collagen Matrix, Inc. and distributed by Biomet 3i, LLC. For additional product information, please refer to the individual product labeling/IFU. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. This material may not be copied or reprinted without the express written consent of Zimmer Biomet. ZB0030 REV A 11/16 ©2016 Zimmer Biomet. All rights reserved.

