

BioMend[®] & BioMend[®] Extend[™]



30 mm x 40 mm

20 mm x 30 mm

15 mm x 20 mm



ABSORBABLE COLLAGEN MEMBRANES

1 Absorbable¹

Eliminates second-stage surgery for membrane removal, designed to reduce wound trauma and surgical chair time

2 Cell - Occlusive

Serves as barrier to prevent epithelial cell migration and allow passage of essential nutrients¹

3 Space - Maintaining

A rigid enough scaffold for tissue regeneration in GTR and GBR procedures



BioMend Membrane being prepared for placement in the sinus cavity



ZIMMER BIOMET
Your progress. Our promise.®

BioMend® & BioMend® Extend™

THE ABSORBABLE MEMBRANE RIGID ENOUGH TO CREATE AND MAINTAIN SPACE

CONVENIENT HANDLING

- Tear resistant, suturable, pliable
- Easy to handle even when hydrated
- Rigid enough to create and maintain space

VERSATILE BARRIER FUNCTION

- Indicated for use in guided tissue and guided bone regeneration procedures
- BioMend Membrane is designed to be absorbed in approximately 8 weeks¹
- BioMend Extend Membrane is designed to be absorbed in approximately 18 weeks¹



Exposed buccal wall defect on #3.



Membrane being prepared for placement.



Membrane cut for placement.



Membrane placed in socket.



Sutures placed to close site.



Final aesthetic result.

Clinical photos ©2013 Gary T. Jones, DMD, Sandhills Oral and Maxillofacial, Fayetteville, NC, USA. All rights reserved. Individual results may vary.

Ordering Information:

Catalog Number	Description	Catalog Number	Description
0103Z	BioMend Membrane, 15 x 20 mm	0140Z	BioMend Extend Membrane, 15 x 20 mm
0105Z	BioMend Membrane, 20 x 30 mm	0141Z	BioMend Extend Membrane, 20 x 30 mm
0107Z	BioMend Membrane, 30 x 40 mm	0142Z	BioMend Extend Membrane, 30 x 40 mm

¹ Data on file with Collagen Matrix, Inc.

Contact us at 1-800-342-5454 or visit zimmerbiometdental.com

Zimmer Biomet Dental
 Global Headquarters
 4555 Riverside Drive
 Palm Beach Gardens, FL 33410
 Phone: +1-561-776-6700
 Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc., and distributed and marketed by Zimmer Biomet Dental (and, in the case of distribution and marketing, its authorized marketing partners). In certain countries, BioMend and BioMend Extend Membranes are manufactured by Collagen Matrix, Inc. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZB0056 REV A 08/17 ©2017 Zimmer Biomet. All rights reserved.

