



Quality & Confidence ALL IN ONE

1 SECURE AND SIMPLE

- Easy and stable screw pick-up
- Reliable driver to screw connection
- Stored and sterilized in one durable kit tray

2 COMPACT

- Autoclavable metal storage tray
- Heavy, stainless steel screw driver handle
- 1.5mmD Micro and 2.0mmD Mini System

3 VERSATILE MODULE – MIX AND MATCH

- Titanium alloy fixation screws in 1.5mmD, 2.0mmD and 2.3mmD
- Color coded drills, drivers and screws
- Choose and add for your individual practice needs at any time

Quality and Confidence – ALL IN ONE

The **Zimmer Screw Fixation System** provides a compact solution for the temporary fixation and stabilization of bone transplants, suitable resorbable and non-resorbable bone replacement materials, and membranes for the alveolar ridge. Two color coded systems in 1.5mmD MICRO (blue) and 2.0mmD MINI (red) make easy and rapid identification of the parts possible and simplifies parts matching. This modular storage system permits individual configuration and the open design aims to ensure convenient access during cleaning and sterilization. Fixation screws and mesh are manufactured from pure titanium or titanium alloy. They are biocompatible, corrosion-proof and nontoxic in the biological environment. They allow imaging virtually free of artifacts.



SECURE & SIMPLE The screws are easily picked up and have a reliable connection to the driver.

Getting Started Made Easy with Our Start-Up Screw Fixation Kit



Start-Up Screw Fixation Kit* (69.01.10Z) features:

- Durable, autoclavable metal storage tray
- Heavy, stainless steel screw driver handle
- Micro System, 1.5mmD
 - Stainless steel screw driver inserts (long & short)
 - Stainless steel pilot drill, 14mmL
 - Stainless steel pilot block drill, 15mmL

SCREW FIXATION KIT – START-UP KIT ORDERING INFORMATION

Contents	Tray	Screw Driver Handle	Screw Driver Insert, Short	Screw Driver Insert, Long	Pilot Drill, Micro, 14mmL	Pilot Block Drill, Micro
Assembled Start-Up Kit*	69.01.10Z					
Components	69.01.11Z	75.23.52Z	75.23.23Z	75.23.19Z	69.01.09Z	69.01.16Z

*Screws are not included in the start-up kit

To learn more about the Zimmer Screw Fixation System, please visit us online at www.zimmerdental.com, or to speak to a sales representative call 1 (800) 854-7019 (in the US) or +1 (760) 929-4300 (outside the US).

This document is intended exclusively for physicians and is not intended for laypersons. Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendation. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice in whole or in part. Please refer to the package inserts for important product information, including, but not limited to, contraindications, warnings, precautions, and adverse effects.

To receive our eNews visit us at http://www.zimmerdental.com/news_eNewsletterSignUp.aspx

For more information about our Products, Regenerative Materials and Educational Opportunities, contact us:

In the U.S. 1 (800) 854-7019
 To fax an order 1 (888) 225-2483
 Outside the U.S. +1 (760) 929-4300
 Australia +61 (0)2 9950 5434 or 1 (800) 241 916
 Canada +1 (905) 567-2073 or 1 (800) 265-0968
 Chile +562 231 5185

China +86 21 22115147
 France +33 (0)1 45 12 35 35
 Germany +49-(0)800-101-64-20
 Israel +972 (0)3 6124242
 Italy +39 0438 37681
 Spain +34 900 800 303



Zimmer Dental
 1900 Aston Avenue
 Carlsbad, CA 92008-7308
 USA

www.zimmerdental.com

©2016 Zimmer Dental Inc. All rights reserved. ZD1587EMEA, Rev. 07/16. Screw Fixation System is manufactured by MEDICON eG and distributed by Zimmer Dental Inc. Please note that not all products are available in every country/region. Please check with a Zimmer Dental representative for availability and additional information.