



## Zimmer® Plastic Temporary Abutments



## soft tissue management with esthetics in mind

### 1 extensive selection

- Available in ready-to-use 1mm or 4mm cuff heights for customization
- Straight and angled options enable selection of the optimal abutment for each individual case
- Conveniently packaged with both processing and retaining screws for use in cement or screw-retained restorations

### 2 uncomplicated and adaptable

- Plastic polymer material that is easily modified
- Chair-side preparation that is both brief and cost-effective

### 3 compatible

- PEEK™ plastic polymer, offered in a natural tooth color, bonds well with common temporary crown materials including acrylic and composite

# preserve and maintain soft tissue contours for better outcomes

Zimmer® Plastic Temporary Abutments are a new, improved and expanded product line for use with the *Tapered Screw-Vent®* and *Screw-Vent®* Implant Systems for 3.5mm, 4.5mm and 5.7mm platforms.



New angled abutments enhance your ability to provide the optimal solution

## benefits of Plastic Temporary Abutments



- A simple, esthetic and long-term solution for soft tissue management during the healing period
- Support for the temporary restoration by preserving and maintaining soft tissue contours to create a natural emergence profile and a foundation for a better esthetic result
- Expanded indication up to 180 days in the mouth

### Plastic Temporary Abutments\*

Part Number	PT341S	PT344S	PT341A	PT344A	PT451S	PT454S	PT451A	PT454A	PT561S	PT564S
Implant Platform	3.5mm	3.5mm	3.5mm	3.5mm	4.5mm	4.5mm	4.5mm	4.5mm	5.7mm	5.7mm
Emergence Profile	4.5mm	4.5mm	4.5mm	4.5mm	5.5mm	5.5mm	5.5mm	5.5mm	6.5mm	6.5mm
Cuff Height	1mm	4mm	1mm	4mm	1mm	4mm	1mm	4mm	1mm	4mm
Angulation	straight	straight	angled	angled	straight	straight	angled	angled	straight	straight

\*Call for availability of product in your area.

To learn more about **Zimmer Plastic Temporary Abutments**, please visit us online at [www.zimmerdental.com](http://www.zimmerdental.com) or to speak to a sales representative, call 1 (800) 854 7019.

To receive our eNews visit us at [http://www.zimmerdental.com/news\\_eNewsLetterSignUp.aspx](http://www.zimmerdental.com/news_eNewsLetterSignUp.aspx)

For more information about our Products, Restorative Materials and Educational Opportunities, contact us:

In the U.S. 1 (800) 854 7019  
 To fax an order 1 (888) 225 2483  
 Outside the U.S. +1 (760) 929 4300  
 Australia +61 (0)2 9950 5434 or 1 (800) 241 916  
 Canada +1 (905) 567 2073 or 1 (800) 265 0968

France +33 (0)1 45 12 35 35  
 Germany +49 (0)761 1 56 47 0  
 Israel +972 (0)3 6124242  
 Italy +39 043 855 5573  
 Spain +34 93 846 05 43



Zimmer Dental  
 1900 Aston Avenue  
 Carlsbad, CA 92008-7308  
 USA

[www.zimmerdental.com](http://www.zimmerdental.com)

©2009 Zimmer Dental Inc. All rights reserved. 12664, Rev. 12/09. Please note that not all products are available in every country/region. Please check with a Zimmer Dental representative for availability and additional information. Plastic Temporary Abutment products are manufactured and distributed by Zimmer Dental Inc.