



Hex-Lock[®] Short Abutment & Restorative System



Posterior restorations made simple, immediate and convenient.

Posterior restorations present many challenges for both clinicians and labs. Limited interocclusal space in which to work, reduced visibility, and of course, time-consuming prep work. But now Zimmer is addressing these challenges with the new *Hex-Lock Short Abutment and Restorative System*, designed to simplify the entire posterior restoration process by being both effective and easy to use. With its shorter cone height, predetermined margins and off-the-shelf convenience, chair and preparation times can be significantly reduced. The restorative protocol is the same as our popular *Hex-Lock Contour Abutment* line, requiring no additional learning curve. And with color-coding to match the emergence profiles, you've got a system that's simple, immediate and convenient.

Hex-Lock Short Abutments feature:

- A shorter profile (4mm cone), specifically designed for the limited interocclusal space of the posterior, and predetermined margins to reduce or eliminate the need for preparations
- Straight, continuous margins around the abutment that allow for easy placement in any orientation
- Abutment collars that are contoured below the margin, allowing space for soft tissue around the implant without compromising strength or esthetics
- An all-inclusive Short Abutment Restorative System that includes all components needed for easy impression-taking and immediate provisionalization, saving valuable chair time and money
- Ideal for use in either one- or two-stage surgical protocols



Compatibility and ease at every level.

Hex-Lock Short Abutments are available for all *Tapered Screw-Vent*® Implants and support all our prosthetic platforms. Available with three emergence profiles, two different cuff heights, a complete restorative system and a host of other features, the new *Hex-Lock* Short Abutment and Restorative System delivers ease at every level.

- Color-coded snap-on impression caps allow for simple abutment-level impressions
- Provisional copings help to create a provisional crown to shape gingival tissue while the final restoration is being fabricated
- Analogs and waxing copings assist the laboratory in fabrication of the final restoration



Hex-Lock Short Restorative System includes healing caps, snap-on impression caps, and analogs, as well as provisional and waxing copings.

Color-Coding for Hex-Lock Short Abutment Impression Caps and Analogs

Color	Emergence Profile
Tan	4.5mmD
Rose	5.5mmD
Yellow	6.5mmD

Hex-Lock Short Abutments

Part Number	SA341	SA342	SA451	SA452	SA561	SA562
Implant Platform	3.5mmD	3.5mmD	4.5mmD	4.5mmD	5.7mmD	5.7mmD
Emergence Profile	4.5mmD	4.5mmD	5.5mmD	5.5mmD	6.5mmD	6.5mmD
Cuff Height	1mm	2mm	1mm	2mm	1mm	2mm

Short Healing Caps

Part Number	SAHEAL4	SAHEAL5	SAHEAL6
Emergence Profile	4.5mmD	5.5mmD	6.5mmD

Short Restorative Kits

(5.5mmD size shown only)

Part Number	SAKIT4	SAKIT5	SAKIT6
Emergence Profile	4.5mmD	5.5mmD	6.5mmD

To learn more about *Hex-Lock* Short Abutments, please visit us online at www.zimmerdental.com, or to speak to a sales representative, call 1 (800) 854-7019.

To receive our eNews visit us at http://www.zimmerdental.com/news_eNewsLetterSignUp.aspx

For more information about our Products and Educational Opportunities, contact us:

In the U.S. 1 (800) 854-7019
 To fax an order 1 (888) 225-2483
 Outside the U.S. +1 (760) 929-4300
 Australia +61 (0)2 9950 5434 or 1 (800) 241 916
 Canada +1 (905) 567-2073 or 1 (800) 265-0968

France +33 (0)1 45 12 35 35
 Germany +49 (0)761 1 56 47 0
 Israel +972 (0)3 6124242
 Italy +39 043 855 5573
 Spain +34 93 846 05 43



Zimmer Dental
 1900 Aston Avenue
 Carlsbad, CA 92008-7308
 USA

www.zimmerdental.com