



# Zfx™ construction components

Product Catalogue 2017

€ 100.–  
New customer  
discount!




# Get connected!

For even more flexibility.

Download now the Zfx library for free!

3shape 



DIRECT SERVER  
CONNECTION

**3shape  
Software Update**  
with latest Zimmer/Zfx  
Implant Libraries compatible with Biomet 3i,  
Camlog, Dentsply Implants, Megagen,  
Nobel Biocare, Straumann, Zimmer ...



**Zfx**

**Zfx**

Download now the Zfx implant library  
for free at: [www.zfx-dental.com](http://www.zfx-dental.com)







# Table of Contents

## Zfx™ construction components

<b>3i-CO Type</b> / Conical® by Biomet 3i .....	7	<b>OS-GS Type</b> / GS & TS System® by Osstem.....	42
<b>3i-CE Type</b> / Certain® by Biomet 3i .....	8	<b>PI-AMP Type</b> / Amplified® by P-I .....	43
<b>3i-IO Type</b> / IOL® by Biomet 3i .....	10	<b>PI-CO Type</b> / Conical Abutment® by P-I.....	44
<b>3i-LO Type</b> / Low Profile® by Biomet 3i .....	11	<b>PI-HE Type</b> / External Hexagon® by P-I .....	45
<b>3i-OS Type</b> / ExHex® by Biomet 3i .....	12	<b>PI-MT Type</b> / Morse Taper® by P-I .....	46
<b>BH-IN Type</b> / Tapered Internal® by BioHorizons.....	14	<b>ST-BO Type</b> / Bone Level® by Straumann Institut .....	48
<b>BR-MU Type</b> / SKY uni.cone® / SKY fast and fixed® by Bredent .....	15	<b>ST-MU Type</b> / Multi-Base® by Straumann Institut .....	50
<b>BR-SK Type</b> / SKY® by Bredent .....	16	<b>ST-OC Type</b> / Octa® (External Octagon) by Straumann Institut .....	51
<b>CA-CA Type</b> / Screw Line® by Camlog.....	18	<b>ST-SY Type</b> / SynOcta® by Straumann Institut .....	52
<b>CA-VA Type</b> / VARIO SR® by Camlog.....	20	<b>SW-GL Type</b> / Global® by Sweden & Martina .....	54
<b>DU-SU Type</b> / Superline & Implantium by Dentium .....	21	<b>SW-OU Type</b> / Outlink <sup>2</sup> ® by Sweden & Martina .....	54
<b>AS-EV Type</b> / EV by Dentsply Implants (Astratech).....	22	<b>SW-PR Type</b> / Premium Kohno® by Sweden & Martina.....	55
<b>AS-OS Type</b> / OsseoSpeed® by Dentsply Implants (Astratech) .....	24	<b>TH-SP Type</b> / SPI® by Thommen Medical .....	56
<b>AS-UN Type</b> / UniAbutment® by Dentsply Implants (Astratech) .....	26	<b>ZI-EZ Type</b> / Eztetic Implant® by Zimmer .....	57
<b>DE-MU Type</b> / Multi-Purpose® by Dentsply Implants .....	27	<b>ZI-SC Type</b> / Tapered Screw Vent® by Zimmer .....	58
<b>DE-XI Type</b> / Frialit-Xive® by Dentsply Implants.....	28	<b>ZI-SH Type</b> / Shoulder Abutment® by Zimmer .....	60
<b>ME-ON Type</b> / AnyOne® by Megagen .....	30	<b>ZI-SW Type</b> / SwissPlus® by Zimmer.....	61
<b>ME-AN Type</b> / AnyRidge® by Megagen.....	30	<b>ZI-TA Type</b> / Tapered Abutment® by Zimmer .....	62
<b>ME-EZ Type</b> / EZ Plus® by Megagen.....	31	<b>Additional component accessories</b> .....	63
<b>ME-RE Type</b> / Rescue Internal® by Megagen .....	32	<b>Instructions for use</b> .....	64
<b>MI-IN Type</b> / Internal Hex® by MIS.....	33	<b>Warranties</b> .....	64
<b>NO-AC Type</b> / Active® by Nobel Biocare .....	34	<b>Technical processing</b> .....	65
<b>NO-BR Type</b> / Brånemark® by Nobel Biocare.....	36	<b>Order form</b> .....	66
<b>NO-MU Type</b> / Multi-unit® by Nobel Biocare.....	38		
<b>NO-RE Type</b> / Replace® by Nobel Biocare.....	40		

# Legend

## Ø Diameter












	3.0 – 3.9 mm
	4.0 – 4.9 mm
	5.0 – 5.9 mm
	6.0 – 6.9 mm

  
**Recommended Torque** Torque for screws

Note: Products which are identified with a "®" are registered trademarks of the corresponding Manufacturer. The representations of the components could differ from the original.

The indications relate to the use of the geometry itself and not exclusively on the titanium base.

## Indications

	Suitable for full metal abutments
	Suitable for full zirconium abutments
	Suitable for zirconium abutments from zirconium for bonding on titanium base
	Suitable for full metal anatomical abutments
	Suitable for full anatomical zirconium
	Suitable for anatomical zirconium abutments for bonding on titanium base
	Suitable for full metal bridges (non-engaging)
	Suitable for full zirconium bridges (non-engaging)
	Suitable for zirconium bridges for bonding on titanium base (non-engaging interface)
	Technically possible, however not recommended. No warranty!
	Not feasible!

## The Zfx goal: To create the ideal emergence profile

Dental implants have become a common restorative option and their use has revolutionized the science of restorative dentistry. The ultimate goal of an implant restoration is to be functional and aesthetic.

The emergence profile plays a crucial role in creating a healthy and natural looking implant supported crown. The final position and shape of the gingival margin can be influenced by changing the contour of the emergence profile and thus modifying the pressure exerted on these tissues.

Optimal aesthetics can be achieved by utilizing the emergence profile to restore the height and shape of the peri-coronal soft tissues whilst long term success is maintained by providing a stable and healthy abutment-to-tissue relationship.

## The Zfx solution: Introducing the Zfx "Small Series" titanium implant interface for hybrid zirconia and CoCr Abutments

Zfx is proud to introduce our new "Small Series" range of titanium interfaces with a reduced implant level diameter, allowing for improved flexibility when designing the emergence profile of implant restorations. The height of the platform on which the final restoration is cemented has also been lowered, creating an even thinner metal margin at the base of the abutment. This has significantly reduced any shadowing showing through the soft tissue, thus enabling the creation of an aesthetically superior abutment.

These "Small Series" range of interfaces are suitable for all indications and provide much needed versatility in the overall design of an implant abutment restoration. They are particularly useful in compromised cases involving thin soft tissue biotypes, tight, fibrous tissue profiles or shallow implant placements where it has been traditionally difficult to achieve ideal results.

The "Regular Series" range of titanium interfaces are still available and are indicated for thicker soft tissue biotypes, where increased soft tissue pressure from a flared emergence profile is indicated to shape and contour the surrounding soft tissue.




Emergence Profile



Regular Design Series



Small Design Series



Regular Design vs. Small Design

## Zfx™ Evolution Matchholder

These Matchholders (scan bodies) have to be used with the laboratory scanners of the Zfx™ Evolution plus. A code system is available on the Matchholders surface; this code is captured during the scan process and automatically detects and recognizes the used implant system and its platform size. Due to its high accuracy, the Matchholder is suitable for scanning implant bridges. All Zfx™ Evolution Matchholders have an integrated screw that have to be utilized with the Torque Wrench for Matchholders. The integrated screws cannot be utilized with the original key of the implant system.



Zfx™ Evolution Matchholder



Torque wrench for Zfx™ Evolution Matchholders, Art.-No.: ZFX02002067

## Torque Wrench for Matchholders

The Torque Wrench is required for tightening the Zfx™ Evolution Matchholders (scan bodies). It is at the same time, with a breadth of 1.25 mm, a universal key. The Torque Wrench also has a slipping clutch that guarantees an equal tightening torque for each match holder. Art-No : ZFX02002067



Zfx™ IntraScan Matchholder

## Zfx™ IntraScan Matchholder


These Matchholders (scan bodies) have to be used with the Zfx™ IntraScan in the oral cavity of patients. The Zfx™ IntraScan Matchholders are available in 2 different heights, H4 (height: 4 mm) or H7 (height: 7 mm). These two heights have to be used at the patient accordingly to the different implant's depth and gingiva's height. The Matchholder H7 has an additional step in order to avoid any confusion between the two heights when designing in the CAD software. All Zfx™ IntraScan Matchholder have an integrated screw, that has to be used with the Hex driver, Long Standard for Zfx™ IntraScan (Art.-No.: HXL1.25) or with the Hex driver, Long with Gemlock Retention for Zfx™ IntraScan (Art.-No.: HXLGR1.25). The integrated screws cannot be utilized with the original key of the implant system.



Hex driver, Long Standard for Zfx™ IntraScan and Hex driver, Long with Gemlock Retention for Zfx™ IntraScan (HXL1.25 & HXLGR1.25)

## Model analogues

The model analogues have to be used for the models, which have been digitally created with the Zfx™ Modelcreator. The model analogues are pinned through the lateral bore: In this way they are fixed in exact alignment and height. The model analogues can also be implemented in a plaster model as usual, using an impression.



Zfx™ IntraScan Matchholder

## Zfx™ Pre-Abutments

With the Zfx™ Pre-Abutment Blanks Zfx provides the appropriate Titanium Blanks (Great 5) with pre-manufactured connecting geometries for all common implant systems. The Pre-Abutment Processing Upgrade for the Zfx™ Inhouse5x offers the possibility to produce individual abutments in a time-saving and cost-effective manner in your own lab.



Model analogs for all current implant systems



Zfx™ Pre-Abutments



## 3i-CO Type

Compatible with the implant system **Conical® by Biomet 3i**.

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
		€ 49.-	€ 19.-	€ 19.-
<b>ø 3.4 mm</b>	-	S-TI-Base-3i-CO-TBUL-3.4	3i-CO-SCTI-3.4-6.0	3i-CO-SCTI-3.4-6.0
Art-No.:	-	ZFX11002086	ZFX09001007	ZFX09001007
<b>ø 4.1 / 5.0 mm</b>	-	S-TI-Base-3i-CO-TBUL-4.1-5.0	3i-CO-SCTI-3.4-6.0	3i-CO-SCTI-3.4-6.0
Art-No.:	-	ZFX11002087	ZFX09001007	ZFX09001007

**TWIN PACK**

### Components Accessories

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
	2 PC. € 139.-	2 PC. € 149.-		2 PC. € 44.-
<b>ø 3.4 mm</b>	EVO-3i-CO-MACH-3.4	INTRA-3i-CO-MACH-3.4-H4	-	MC-Analog-3i-CO-3.4
Art-No.:	ZFX05002158	ZFX05001925	-	ZFX15002289
<b>ø 4.1 / 5.0 mm</b>	EVO-3i-CO-MACH-4.1-5.0	INTRA-3i-CO-MACH-4.1-5.0-H4	-	MC-Analog-3i-CO-4.1-5.0
Art-No.:	ZFX05002159	ZFX05001926	-	ZFX15002290

## 3i-CE Type

Compatible with the implant system **Certain®** by **Biomet 3i**.

### Single Implant Restorations

	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
	 € 49.-	 € 49.-	 € 19.-	 € 19.-
<b>ø 3.4 mm</b>	TI-Base-3i-CE-TBLO-3.4	S-TI-Base-3i-CE-TBLO-3.4	3i-CE-SCTI-3.4-6.0	3i-CE-SCZI-3.4-6.0
Art-No.:	ZFX11000668	ZFX11001096	ZFX09000978	ZFX09000610
<b>ø 4.1 mm</b>	TI-Base-3i-CE-TBLO-4.1	S-TI-Base-3i-CE-TBLO-4.1	3i-CE-SCTI-3.4-6.0	3i-CE-SCZI-3.4-6.0
Art-No.:	ZFX11000669	ZFX11001097	ZFX09000978	ZFX09000610
<b>ø 5.0 mm</b>	TI-Base-3i-CE-TBLO-5.0	S-TI-Base-3i-CE-TBLO-5.0	3i-CE-SCTI-3.4-6.0	3i-CE-SCZI-3.4-6.0
Art-No.:	ZFX11000670	ZFX11001098	ZFX09000978	ZFX09000610
<b>ø 6.0 mm</b>	-	S-TI-Base-3i-CE-TBLO-6.0	3i-CE-SCTI-3.4-6.0	3i-CE-SCZI-3.4-6.0
Art-No.:	-	ZFX11000671	ZFX09000978	ZFX09000610





### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
	 € 49.-	 € 49.-	 € 19.-	 € 19.-
<b>ø 3.4 mm</b>	TI-Base-3i-CE-TBUL-3.4	S-TI-Base-3i-CE-TBUL-3.4	3i-CE-SCTI-3.4-6.0	3i-CE-SCZI-3.4-6.0
Art-No.:	ZFX11000672	ZFX11001100	ZFX09000978	ZFX09000610
<b>ø 4.1 mm</b>	TI-Base-3i-CE-TBUL-4.1	S-TI-Base-3i-CE-TBUL-4.1	3i-CE-SCTI-3.4-6.0	3i-CE-SCZI-3.4-6.0
Art-No.:	ZFX11000673	ZFX11001101	ZFX09000978	ZFX09000610
<b>ø 5.0 mm</b>	TI-Base-3i-CE-TBUL-5.0	S-TI-Base-3i-CE-TBUL-5.0	3i-CE-SCTI-3.4-6.0	3i-CE-SCZI-3.4-6.0
Art-No.:	ZFX11000674	ZFX11001102	ZFX09000978	ZFX09000610
<b>ø 6.0 mm</b>	-	S-TI-Base-3i-CE-TBUL-6.0	3i-CE-SCTI-3.4-6.0	3i-CE-SCZI-3.4-6.0
Art-No.:	-	ZFX11000675	ZFX09000978	ZFX09000610



TWIN PACK

Components Accessories

	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.4 mm</b>	EVO-3I-CE-MACH-3.4	INTRA-3I-CE-MACH-3.4-H4	INTRA-3I-CE-MACH-3.4-H7	MC-Analog-3I-CE-3.4
Art-No.:	ZFX05002018	ZFX05001942	ZFX05001943	ZFX15002116
<b>ø 4.1 mm</b>	EVO-3I-CE-MACH-4.1	INTRA-3I-CE-MACH-4.1-H4	INTRA-3I-CE-MACH-4.1-H7	MC-Analog-3I-CE-4.1
Art-No.:	ZFX05002019	ZFX05001944	ZFX05001945	ZFX15002117
<b>ø 5.0 mm</b>	EVO-3I-CE-MACH-5.0	INTRA-3I-CE-MACH-5.0-H4	INTRA-3I-CE-MACH-5.0-H7	MC-Analog-3I-CE-5.0
Art-No.:	ZFX05002020	ZFX05001946	ZFX05001947	ZFX15002118
<b>ø 6.0 mm</b>	EVO-3I-CE-MACH-6.0	INTRA-3I-CE-MACH-6.0-H4	INTRA-3I-CE-MACH-6.0-H7	MC-Analog-3I-CE-6.0
Art-No.:	ZFX05002021	ZFX05001948	ZFX05001949	ZFX15002119

Zfx™ Pre-Abutment Blanks including screw



<b>ø 3.4 mm</b>	Art-No.:	ZFX18000350
<b>ø 4.1 - 6.0 mm</b>	Art-No.:	ZFX18000360

## 3i-IO Type

Compatible with the implant system **IOL® by Biomet 3i**.

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
		€ 49.-	€ 19.-	€ 19.-
<b>ø 4.5 mm</b>	-	S-TI-Base-3i-IO-TBUL-4.5	3i-CO-SCTI-3.4-6.0	3i-CO-SCTI-3.4-6.0
Art-No.:	-	ZFX11002089	ZFX09001007	ZFX09001007
			10 Ncm	10 Ncm

TWIN PACK

### Components Accessories

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
	2 PC. € 139.-	2 PC. € 149.-		2 PC. € 44.-
<b>ø 4.5 mm</b>	EVO-3i-IO-MACH-4.5	INTRA-3i-IO-MACH-4.5-H4	-	MC-Analog-3i-IO-4.5
Art-No.:	ZFX05002161	ZFX05001927	-	ZFX15002291

# 3i-LO Type

Compatible with the implant system **Low Profile®** by Biomet 3i.

## Implant Bar / Bridge Restorations

	<i>Titanium base not protected against rotation - Regular Design</i>	<i>Titanium base not protected against rotation - Small Design</i>	<i>Bridge screw Standard</i>	<i>Bridge screw for All-Zirconia Implant Bridges</i>
		€ 49.-	€ 19.-	€ 19.-
<b>ø 4.8 mm</b>	-	S-TI-Base-3I-LO-TBUL-4.8	3I-LO-SCTI-4.8	3I-LO-SCTI-4.8
Art-No.:	-	ZFX11002090	ZFX09001008	ZFX09001008
			10 Ncm	10 Ncm

TWIN PACK

## Components Accessories

	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	2 PC. € 139.-	2 PC. € 149.-		2 PC. € 44.-
<b>ø 4.8 mm</b>	EVO-3I-LO-MACH-4.8	INTRA-3I-LO-MACH-4.8-H4	-	MC-Analog-3I-LO-4.8
Art-No.:	ZFX05002162	ZFX05001928	-	ZFX15002292





## 3i-OS Type

Compatible with the implant system **ExHex® by Biomet 3i**.

### Single Implant Restorations

	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
	 € 49.-	 € 49.-	 € 19.-	 € 19.-
<b>ø 3.4 mm</b>	TI-Base-3i-OS-TBLO-3.4	S-TI-Base-3i-OS-TBLO-3.4	3i-OS-SCTI-3.4-6.0	3i-OS-SCZI-3.4-6.0
Art-No.:	ZFX11000676	ZFX11001104	ZFX09000979	ZFX09000611
<b>ø 4.1 mm</b>	TI-Base-3i-OS-TBLO-4.1	S-TI-Base-3i-OS-TBLO-4.1	3i-OS-SCTI-3.4-6.0	3i-OS-SCZI-3.4-6.0
Art-No.:	ZFX11000677	ZFX11001105	ZFX09000979	ZFX09000611
<b>ø 5.0 mm</b>	TI-Base-3i-OS-TBLO-5.0	S-TI-Base-3i-OS-TBLO-5.0	3i-OS-SCTI-3.4-6.0	3i-OS-SCZI-3.4-6.0
Art-No.:	ZFX11000678	ZFX11001106	ZFX09000979	ZFX09000611
<b>ø 6.0 mm</b>	-	S-TI-Base-3i-OS-TBLO-6.0	3i-OS-SCTI-3.4-6.0	3i-OS-SCZI-3.4-6.0
Art-No.:	-	ZFX11000679	ZFX09000979	ZFX09000611

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
		 € 49.-	 € 19.-	 € 19.-
<b>ø 3.4 mm</b>	-	S-TI-Base-3i-OS-TBUL-3.4	3i-OS-SCTI-3.4-6.0	3i-OS-SCZI-3.4-6.0
Art-No.:	-	ZFX11000680	ZFX09000979	ZFX09000611
<b>ø 4.1 mm</b>	-	S-TI-Base-3i-OS-TBUL-4.1	3i-OS-SCTI-3.4-6.0	3i-OS-SCZI-3.4-6.0
Art-No.:	-	ZFX11000681	ZFX09000979	ZFX09000611
<b>ø 5.0 mm</b>	-	S-TI-Base-3i-OS-TBUL-5.0	3i-OS-SCTI-3.4-6.0	3i-OS-SCZI-3.4-6.0
Art-No.:	-	ZFX11000682	ZFX09000979	ZFX09000611
<b>ø 6.0 mm</b>	-	S-TI-Base-3i-OS-TBUL-6.0	3i-OS-SCTI-3.4-6.0	3i-OS-SCZI-3.4-6.0
Art-No.:	-	ZFX11000683	ZFX09000979	ZFX09000611

TWIN PACK

Components Accessories

	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.4 mm</b>	EVO-3I-OS-MACH-3.4	INTRA-3I-OS-MACH-3.4-H4	INTRA-3I-OS-MACH-3.4-H7	MC-Analog-3I-OS-3.4
Art-No.:	ZFX05002022	ZFX05001950	ZFX05001951	ZFX15002120
<b>ø 4.1 mm</b>	EVO-3I-OS-MACH-4.1	INTRA-3I-OS-MACH-4.1-H4	INTRA-3I-OS-MACH-4.1-H7	MC-Analog-3I-OS-4.1
Art-No.:	ZFX05002023	ZFX05001952	ZFX05001953	ZFX15002121
<b>ø 5.0 mm</b>	EVO-3I-OS-MACH-5.0	INTRA-3I-OS-MACH-5.0-H4	INTRA-3I-OS-MACH-5.0-H7	MC-Analog-3I-OS-5.0
Art-No.:	ZFX05002024	ZFX05001954	ZFX05001955	ZFX15002122
<b>ø 6.0 mm</b>	EVO-3I-OS-MACH-6.0	INTRA-3I-OS-MACH-6.0-H4	INTRA-3I-OS-MACH-6.0-H7	MC-Analog-3I-OS-6.0
Art-No.:	ZFX05002025	ZFX05001956	ZFX05001957	ZFX15002123

Zfx™ Pre-Abutment Blanks including screw







<b>ø 3.4 mm</b>	Art-No.:	ZFX18000370
<b>ø 4.1 - 6.0 mm</b>	Art-No.:	ZFX18000380

## BH-IN Type

Compatible with the implant system **Tapered Internal®** by BioHorizons.

**TWIN PACK**

### Components Accessories






	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.5 mm</b>	EVO-BH-IN-MACH-3.5	INTRA-BH-IN-MACH-3.5-H4	INTRA-BH-IN-MACH-3.5-H7	MC-Analog-BH-IN-3.5
Art-No.:	ZFX05002282	ZFX05002285	ZFX05002288	ZFX15002331
<b>ø 4.5 mm</b>	EVO-BH-IN-MACH-4.5	INTRA-BH-IN-MACH-4.5-H4	INTRA-BH-IN-MACH-4.5-H7	MC-Analog-BH-IN-4.5
Art-No.:	ZFX05002283	ZFX05002286	ZFX05002289	ZFX15002332
<b>ø 5.7 mm</b>	EVO-BH-IN-MACH-5.7	INTRA-BH-IN-MACH-5.7-H4	INTRA-BH-IN-MACH-5.7-H7	MC-Analog-BH-IN-5.7
Art-No.:	ZFX05002284	ZFX05002287	ZFX05002290	ZFX15002333



## BR-MU Type





Compatible with the implant system **SKY uni.cone® / SKY fast and fixed®** by Bredent.

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
		 € 49.-	 € 19.-	 € 19.-
<b>ø 4.5 mm</b>	-	S-TI-Base-BR-MU-TBUL-4.5	BR-MU-SCTI-4.5-5.6	BR-MU-SCTI-4.5-5.6
Art-No.:	-	ZFX11002199	ZFX09002186	ZFX09002186
<b>ø 5.6 mm</b>	-	S-TI-Base-BR-MU-TBUL-5.6	BR-MU-SCTI-4.5-5.6	BR-MU-SCTI-4.5-5.6
Art-No.:	-	ZFX11002196	ZFX09002186	ZFX09002186

TWIN PACK






### Components Accessories

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
 2 PC. € 139.-	 2 PC. € 149.-		 2 PC. € 44.-	
<b>ø 4.5 mm</b>	EVO-BR-MU-MACH-4.5	INTRA-BR-MU-MACH-4.5-H4	-	MC-Analog-BR-MU-4.5
Art-No.:	ZFX05002061	ZFX05001150	-	ZFX15002329
<b>ø 5.6 mm</b>	EVO-BR-MU-MACH-5.6	INTRA-BR-MU-MACH-5.6-H4	-	MC-Analog-BR-MU-5.6
Art-No.:	ZFX05002179	ZFX05001939	-	ZFX15002330






## BR-SK Type

Compatible with the implant system **SKY® by Bredent.**

### Single Implant Restorations





	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
		 € 49.-	 € 19.-	 € 19.-
<b>ø 3.5 mm</b>	-	S-TI-Base-BR-SK-TBLO-3.5	BR-SK-SCTI-3.5-4.0	BR-SK-SCTI-3.5-4.0
Art-No.:	-	ZFX11002205	ZFX09002185	ZFX09002185
<b>ø 4.0 mm</b>	-	S-TI-Base-BR-SK-TBLO-4.0	BR-SK-SCTI-3.5-4.0	BR-SK-SCTI-3.5-4.0
Art-No.:	-	ZFX11002194	ZFX09002185	ZFX09002185
			↑ Recommended Torque	↑ Recommended Torque
			25 Ncm	25 Ncm

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
		 € 49.-	 € 19.-	 € 19.-
<b>ø 3.5 mm</b>	-	S-TI-Base-BR-SK-TBUL-3.5	BR-SK-SCTI-3.5-4.0	BR-SK-SCTI-3.5-4.0
Art-No.:	-	ZFX11002207	ZFX09002185	ZFX09002185
<b>ø 4.0 mm</b>	-	S-TI-Base-BR-SK-TBUL-4.0	BR-SK-SCTI-3.5-4.0	BR-SK-SCTI-3.5-4.0
Art-No.:	-	ZFX11002195	ZFX09002185	ZFX09002185
			↑ Recommended Torque	↑ Recommended Torque
			25 Ncm	25 Ncm

TWIN PACK

Components Accessories

	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.5 mm</b>	EVO-BR-SK-MACH-3.5	INTRA-BR-SK-MACH-3.5-H4	INTRA-BR-SK-MACH-3.5-H7	MC-Analog-BR-SK-3.5
Art-No.:	ZFX05002256	ZFX05001132	ZFX05001133	ZFX15002327
<b>ø 4.0 mm</b>	EVO-BR-SK-MACH-4.0	INTRA-BR-SK-MACH-4.0-H4	INTRA-BR-SK-MACH-4.0-H7	MC-Analog-BR-SK-4.0
Art-No.:	ZFX05002178	ZFX05001102	ZFX05001103	ZFX15002328

## CA-CA Type

Compatible with the implant system **Screw Line®** by Camlog.

### Single Implant Restorations

	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
			<b>T</b> Recommended Torque	<b>T</b> Recommended Torque
<b>ø 3.3 mm</b>	-	-	CA-CA-SCTI-3.3-4.3	CA-CA-SCTI-3.3-4.3
Art-No.:	-	K2244.3348*	ZFX09002323	ZFX09002323
<b>ø 3.8 mm</b>	-	-	CA-CA-SCTI-3.3-4.3	CA-CA-SCTI-3.3-4.3
Art-No.:	-	K2244.3848*	ZFX09002323	ZFX09002323
<b>ø 4.3 mm</b>	-	-	CA-CA-SCTI-3.3-4.3	CA-CA-SCTI-3.3-4.3
Art-No.:	-	K2244.4348*	ZFX09002323	ZFX09002323
<b>ø 5.0 mm</b>	-	-	CA-CA-SCTI-5.0-6.0	CA-CA-SCTI-5.0-6.0
Art-No.:	-	K2244.5048*	ZFX09002324	ZFX09002324
<b>ø 6.0 mm</b>	-	-	CA-CA-SCTI-5.0-6.0	CA-CA-SCTI-5.0-6.0
Art-No.:	-	K2244.6048*	ZFX09002324	ZFX09002324





### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
			<b>T</b> Recommended Torque	<b>T</b> Recommended Torque
<b>ø 3.3 mm</b>	-	S-TI-Base-CA-CA-TBUL-3.3	CA-CA-SCTI-3.3-4.3	CA-CA-SCTI-3.3-4.3
Art-No.:	-	ZFX11002208	ZFX09002323	ZFX09002323
<b>ø 3.8 mm</b>	-	S-TI-Base-CA-CA-TBUL-3.8	CA-CA-SCTI-3.3-4.3	CA-CA-SCTI-3.3-4.3
Art-No.:	-	ZFX11002209	ZFX09002323	ZFX09002323
<b>ø 4.3 mm</b>	-	S-TI-Base-CA-CA-TBUL-4.3	CA-CA-SCTI-3.3-4.3	CA-CA-SCTI-3.3-4.3
Art-No.:	-	ZFX11002210	ZFX09002323	ZFX09002323
<b>ø 5.0 mm</b>	-	S-TI-Base-CA-CA-TBUL-5.0	CA-CA-SCTI-5.0-6.0	CA-CA-SCTI-5.0-6.0
Art-No.:	-	ZFX11002211	ZFX09002324	ZFX09002324
<b>ø 6.0 mm</b>	-	S-TI-Base-CA-CA-TBUL-6.0	CA-CA-SCTI-5.0-6.0	CA-CA-SCTI-5.0-6.0
Art-No.:	-	ZFX11002212	ZFX09002324	ZFX09002324

\*Not available at Zfx. This code is the original Camlog Institute's article code.

TWIN PACK

Components Accessories

	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.3 mm</b>	EVO-CA-CA-MACH-3.3	INTRA-CA-CA-MACH-3.3-H4	INTRA-CA-CA-MACH-3.3-H7	MC-Analog-CA-CA-3.3
Art-No.:	ZFX05002029	ZFX05001964	ZFX05001965	ZFX15002127
<b>ø 3.8 mm</b>	EVO-CA-CA-MACH-3.8	INTRA-CA-CA-MACH-3.8-H4	INTRA-CA-CA-MACH-3.8-H7	MC-Analog-CA-CA-3.8
Art-No.:	ZFX05002030	ZFX05001966	ZFX05001967	ZFX15002128
<b>ø 4.3 mm</b>	EVO-CA-CA-MACH-4.3	INTRA-CA-CA-MACH-4.3-H4	INTRA-CA-CA-MACH-4.3-H7	MC-Analog-CA-CA-4.3
Art-No.:	ZFX05002031	ZFX05001968	ZFX05001969	ZFX15002129
<b>ø 5.0 mm</b>	EVO-CA-CA-MACH-5.0	INTRA-CA-CA-MACH-5.0-H4	INTRA-CA-CA-MACH-5.0-H7	MC-Analog-CA-CA-5.0
Art-No.:	ZFX05002032	ZFX05001970	ZFX05001971	ZFX15002130
<b>ø 6.0 mm</b>	EVO-CA-CA-MACH-6.0	INTRA-CA-CA-MACH-6.0-H4	INTRA-CA-CA-MACH-6.0-H7	MC-Analog-CA-CA-6.0
Art-No.:	ZFX05002033	ZFX05001972	ZFX05001973	ZFX15002131

## CA-VA Type

Compatible with the implant system **VARIO SR®** by Camlog.

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
		€ 49.-	€ 19.-	€ 19.-
<b>ø 3.8 / 4.3 mm</b>	-	S-TI-Base-CA-VA-TBUL-3.8-4.3	CA-VA-SCTI-3.8-6.0	15 Ncm
Art-No.:	-	ZFX11002197	ZFX09002187	ZFX09002187
<b>ø 5.0 / 6.0 mm</b>	-	S-TI-Base-CA-VA-TBUL-5.0-6.0	CA-VA-SCTI-3.8-6.0	15 Ncm
Art-No.:	-	ZFX11002198	ZFX09002187	ZFX09002187


### Components Accessories TWIN PACK

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
	2 PC. € 139.-	2 PC. € 149.-		2 PC. € 44.-
<b>ø 3.8 / 4.3 mm</b>	EVO-CA-VA-MACH-3.8-4.3	INTRA-CA-VA-MACH-3.8-4.3-H4	-	MC-Analog-CA-VA-3.8-4.3
Art-No.:	ZFX05002180	ZFX05001940	-	ZFX15002295
<b>ø 5.0 / 6.0 mm</b>	EVO-CA-VA-MACH-5.0-6.0	INTRA-CA-VA-MACH-5.0-6.0-H4	-	MC-Analog-CA-VA-5.0-6.0
Art-No.:	ZFX05002181	ZFX05001941	-	ZFX15002296



## DU-SU Type






Compatible with the implant system **Superline & Implantium®** by Dentium.

Components Accessories		<i>TWIN PACK</i>			
		<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
 <p>2 PC. € 139.-</p>					
<b>ø 3.6/7.0 mm</b>	EVO-DU-SU-MACH-3.6-7.0	-	-	-	
Art-No.:	ZFX05002200	-	-	-	






## AS-EV Type

Compatible with the implant system **EV by Dentsply Implants (Astratech)**.

### Single Implant Restorations


	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
			 € 19.–	 € 19.–
			<b>Recommended Torque</b>	<b>Recommended Torque</b>
<b>ø 3.0 mm</b>	-	-	AS-EV-SCTI-3.0	AS-OS-SCTI-3.0
Art-No.:	-	-	ZFX09002330	ZFX09002330
<b>ø 3.6 mm</b>	-	-	AS-EV-SCTI-3.6	AS-EV-SCTI-3.6
Art-No.:	-	-	ZFX09002331	ZFX09002331
<b>ø 4.2 mm</b>	-	-	AS-EV-SCTI-4.2	AS-EV-SCTI-4.2
Art-No.:	-	-	ZFX09002332	ZFX09002332
<b>ø 4.8 / 5.4 mm</b>	-	-	AS-EV-SCTI-4.8-5.4	AS-EV-SCTI-4.8-5.4
Art-No.:	-	-	ZFX09002333	ZFX09002333

### Implant Bar / Bridge Restorations


	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
			 € 19.–	 € 19.–
			<b>Recommended Torque</b>	<b>Recommended Torque</b>
<b>ø 3.0 mm</b>	-	-	AS-EV-SCTI-3.0	AS-OS-SCTI-3.0
Art-No.:	-	-	ZFX09002330	ZFX09002330
<b>ø 3.6 mm</b>	-	-	AS-EV-SCTI-3.6	AS-EV-SCTI-3.6
Art-No.:	-	-	ZFX09002331	ZFX09002331
<b>ø 4.2 mm</b>	-	-	AS-EV-SCTI-4.2	AS-EV-SCTI-4.2
Art-No.:	-	-	ZFX09002332	ZFX09002332
<b>ø 4.8 / 5.4 mm</b>	-	-	AS-EV-SCTI-4.8-5.4	AS-EV-SCTI-4.8-5.4
Art-No.:	-	-	ZFX09002333	ZFX09002333

Components Accessories

TWIN PACK

	Evolution Matchholder	IntraScan Matchholder H4	IntraScan Matchholder H7	Model analogue
 <p>2 PC. € 139.-</p>				
<b>ø 3.0 mm</b>	EVO-AS-EV-MACH-3.0	-	-	-
Art-No.:	ZFX05002186	-	-	-
<b>ø 3.6 mm</b>	EVO-AS-EV-MACH-3.6	-	-	-
Art-No.:	ZFX05002188	-	-	-
<b>ø 4.2 mm</b>	EVO-AS-EV-MACH-4.2	-	-	-
Art-No.:	ZFX05002190	-	-	-
<b>ø 4.8 mm</b>	EVO-AS-EV-MACH-4.8	-	-	-
Art-No.:	ZFX05002192	-	-	-
<b>ø 5.4 mm</b>	EVO-AS-EV-MACH-5.4	-	-	-
Art-No.:	ZFX05002194	-	-	-








Zfx™ Pre-Abutment Blanks including screw

 <p>€ 44.90</p>	<b>ø 3.0 mm</b>	Art-No.:	ZFX18000300
	<b>ø 3.6 mm</b>	Art-No.:	ZFX18000310
	<b>ø 4.2 mm</b>	Art-No.:	ZFX18000320
	<b>ø 4.8 mm</b>	Art-No.:	ZFX18000330
	<b>ø 5.4 mm</b>	Art-No.:	ZFX18000340








## AS-OS Type

Compatible with the implant system **Osseospeed® by Dentsply Implants (Astratech)**.

### Single Implant Restorations





	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
				
	€ 49.-	€ 49.-	€ 19.-	€ 19.-
			 Recommended Torque	 Recommended Torque
<b>ø 3.0 mm</b>	-	S-TI-Base-AS-OS-TBLO-3.0	AS-OS-SCTI-3.0	AS-OS-SCTI-3.0
Art-No.:	-	ZFX11000684	ZFX09000612	ZFX09000612
<b>ø 3.5 / 4.0 mm</b>	TI-Base-AS-OS-TBLO-3.5-4.0	S-TI-Base-AS-OS-TBLO-3.5-4.0	AS-OS-SCTI-3.5-4.0	AS-OS-SCTI-3.5-4.0
Art-No.:	ZFX11000685	ZFX11001113	ZFX09000613	ZFX09000613
<b>ø 4.5 / 5.0 mm</b>	TI-Base-AS-OS-TBLO-4.5-5.0	S-TI-Base-AS-OS-TBLO-4.5-5.0	AS-OS-SCTI-4.5-5.0	AS-OS-SCTI-4.5-5.0
Art-No.:	ZFX11000686	ZFX11001114	ZFX09000614	ZFX09000614

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
				
	€ 49.-	€ 49.-	€ 19.-	€ 19.-
			 Recommended Torque	 Recommended Torque
<b>ø 3.0 mm</b>	-	S-TI-Base-AS-OS-TBUL-3.0	AS-OS-SCTI-3.0	AS-OS-SCTI-3.0
Art-No.:	-	ZFX11000687	ZFX09000612	ZFX09000612
<b>ø 3.5 / 4.0 mm</b>	TI-Base-AS-OS-TBUL-3.5-4.0	S-TI-Base-AS-OS-TBUL-3.5-4.0	AS-OS-SCTI-3.5-4.0	AS-OS-SCTI-3.5-4.0
Art-No.:	ZFX11000688	ZFX11001116	ZFX09000613	ZFX09000613
<b>ø 4.5 / 5.0 mm</b>	TI-Base-AS-OS-TBUL-4.5-5.0	S-TI-Base-AS-OS-TBUL-4.5-5.0	AS-OS-SCTI-4.5-5.0	AS-OS-SCTI-4.5-5.0
Art-No.:	ZFX11000689	ZFX11001117	ZFX09000614	ZFX09000614

TWIN PACK

Components Accessories

	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.0 mm</b>	EVO-AS-OS-MACH-3.0	INTRA-AS-OS-MACH-3.0-H4	INTRA-AS-OS-MACH-3.0-H7	MC-Analog-AS-OS-3.0
Art-No.:	ZFX05002026	ZFX05001958	ZFX05001959	ZFX15002124
<b>ø 3.5 / 4.0 mm</b>	EVO-AS-OS-MACH-3.5-4.0	INTRA-AS-OS-MACH-3.5-4.0-H4	INTRA-AS-OS-MACH-3.5-4.0-H7	MC-Analog-AS-OS-3.5-4.0
Art-No.:	ZFX05002027	ZFX05001960	ZFX05001961	ZFX15002125
<b>ø 4.5 / 5.0 mm</b>	EVO-AS-OS-MACH-4.5-5.0	INTRA-AS-OS-MACH-4.5-5.0-H4	INTRA-AS-OS-MACH-4.5-5.0-H7	MC-Analog-AS-OS-4.5-5.0
Art-No.:	ZFX05002028	ZFX05001962	ZFX05001963	ZFX15002126

Zfx™ Pre-Abutment Blanks including screw



<b>ø 3.0 mm</b>	Art-No.:	ZFX18000140
<b>ø 3.5 - 4.0 mm</b>	Art-No.:	ZFX18000150
<b>ø 4.5 - 5.0 mm</b>	Art-No.:	ZFX18000160

## AS-UN Type

Compatible with the implant system **UniAbutment®** by Dentsply Implants (Astratech).

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
		€ 49.-	€ 19.-	€ 19.-
<b>20°</b>	-	S-TI-Base-AS-UN-TBUL-20°	AS-UN-SCTI-20°-45°	AS-UN-SCTI-20°-45°
Art-No.:	-	ZFX11002091	ZFX09001009	ZFX09001009
<b>45°</b>	-	S-TI-Base-AS-UN-TBUL-45°	AS-UN-SCTI-20°-45°	AS-UN-SCTI-20°-45°
Art-No.:	-	ZFX11002092	ZFX09001009	ZFX09001009

### Components Accessories TWIN PACK

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
	2 PC. € 139.-	2 PC. € 149.-		2 PC. € 44.-
<b>20°</b>	EVO-AS-UN-MACH-20°	INTRA-AS-UN-MACH-20°-H4	-	MC-Analog-AS-UN-20°
Art-No.:	ZFX05002163	ZFX05001929	-	ZFX15002293
<b>45°</b>	EVO-AS-UN-MACH-45°	INTRA-AS-UN-MACH-45°-H4	-	MC-Analog-AS-UN-45°
Art-No.:	ZFX05002164	ZFX05001930	-	ZFX15002294



# DE-MU Type

Compatible with the implant system **Multi-Purpose®** by Dentsply Implants.

## Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
		€ 49.-	€ 19.-	€ 19.-
<b>ø 3.4 mm</b>	-	S-TI-Base-DE-MU-TBUL-3.4-3.8	DE-MU-SCTI-3.4-5.5	DE-MU-SCTI-3.4-5.5
Art-No.:	-	ZFX11002093	ZFX09001010	ZFX09001010
<b>ø 4.5 mm</b>	-	S-TI-Base-DE-MU-TBUL-4.5	DE-MU-SCTI-3.4-5.5	DE-MU-SCTI-3.4-5.5
Art-No.:	-	ZFX11002094	ZFX09001010	ZFX09001010
<b>ø 5.5 mm</b>	-	S-TI-Base-DE-MU-TBUL-5.5	DE-MU-SCTI-3.4-5.5	DE-MU-SCTI-3.4-5.5
Art-No.:	-	ZFX11002095	ZFX09001010	ZFX09001010

**TWIN PACK**






## Components Accessories

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
	2 PC. € 139.-	2 PC. € 149.-		2 PC. € 44.-
<b>ø 3.4 mm</b>	EVO-DE-MU-MACH-3.4-3.8	INTRA-DE-MU-MACH-3.4-3.8-H4	-	MC-Analog-DE-MU-3.4-3.8
Art-No.:	ZFX05002165	ZFX05001931	-	ZFX15002297
<b>ø 4.5 mm</b>	EVO-DE-MU-MACH-4.5	INTRA-DE-MU-MACH-4.5-H4	-	MC-Analog-DE-MU-4.5
Art-No.:	ZFX05002166	ZFX05001932	-	ZFX15002298
<b>ø 5.5 mm</b>	EVO-DE-MU-MACH-5.5	INTRA-DE-MU-MACH-5.5-H4	-	MC-Analog-DE-MU-5.5
Art-No.:	ZFX05002167	ZFX05001933	-	ZFX15002299






## DE-XI Type

Compatible with the implant system **Frialit-Xive®** by **Dentsply Implants**.

### Single Implant Restorations





	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
				
	€ 49.-	€ 49.-	€ 19.-	€ 19.-
			<b>Recommended Torque</b>	<b>Recommended Torque</b>
<b>ø 3.0 mm</b>	-	S-TI-Base-DE-XI-TBLO-3.0	DE-XI-SCTI-3.0	DE-XI-SCZI-3.0
Art-No.:	-	ZFX11000696	ZFX09000620	ZFX09000622
<b>ø 3.4 mm</b>	TI-Base-DE-XI-TBLO-3.4	S-TI-Base-DE-XI-TBLO-3.4	DE-XI-SCTI-3.4-6.5	DE-XI-SCZI-3.4-6.5
Art-No.:	ZFX11000697	ZFX11001119	ZFX09000621	ZFX09000623
<b>ø 3.8 mm</b>	TI-Base-DE-XI-TBLO-3.8	S-TI-Base-DE-XI-TBLO-3.8	DE-XI-SCTI-3.4-6.5	DE-XI-SCZI-3.4-6.5
Art-No.:	ZFX11000698	ZFX11001120	ZFX09000621	ZFX09000623
<b>ø 4.5 mm</b>	TI-Base-DE-XI-TBLO-4.5	S-TI-Base-DE-XI-TBLO-4.5	DE-XI-SCTI-3.4-6.5	DE-XI-SCZI-3.4-6.5
Art-No.:	ZFX11000699	ZFX11001121	ZFX09000621	ZFX09000623
<b>ø 5.5 mm</b>	TI-Base-DE-XI-TBLO-5.5	S-TI-Base-DE-XI-TBLO-5.5	DE-XI-SCTI-3.4-6.5	DE-XI-SCZI-3.4-6.5
Art-No.:	ZFX11000700	ZFX11001122	ZFX09000621	ZFX09000623
<b>ø 6.5 mm</b>	-	S-TI-Base-DE-XI-TBLO-6.5	DE-XI-SCTI-3.4-6.5	DE-XI-SCZI-3.4-6.5
Art-No.:	-	ZFX11000701	ZFX09000621	ZFX09000623

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
				
	€ 49.-	€ 49.-	€ 19.-	€ 19.-
			<b>Recommended Torque</b>	<b>Recommended Torque</b>
<b>ø 3.0 mm</b>	-	S-TI-Base-DE-XI-TBUL-3.0	DE-XI-SCTIUN-3.0	DE-XI-SCZI-3.0
Art-No.:	-	ZFX11000702	ZFX09000980	ZFX09000622
<b>ø 3.4 mm</b>	TI-Base-DE-XI-TBUL-3.4	S-TI-Base-DE-XI-TBUL-3.4	DE-XI-SCTIUN-3.4-6.5	DE-XI-SCZI-3.4-6.5
Art-No.:	ZFX11000703	ZFX11001125	ZFX09000981	ZFX09000623
<b>ø 3.8 mm</b>	TI-Base-DE-XI-TBUL-3.8	S-TI-Base-DE-XI-TBUL-3.8	DE-XI-SCTIUN-3.4-6.5	DE-XI-SCZI-3.4-6.5
Art-No.:	ZFX11000704	ZFX11001126	ZFX09000981	ZFX09000623
<b>ø 4.5 mm</b>	TI-Base-DE-XI-TBUL-4.5	S-TI-Base-DE-XI-TBUL-4.5	DE-XI-SCTIUN-3.4-6.5	DE-XI-SCZI-3.4-6.5
Art-No.:	ZFX11000705	ZFX11001127	ZFX09000981	ZFX09000623
<b>ø 5.5 mm</b>	TI-Base-DE-XI-TBUL-5.5	S-TI-Base-DE-XI-TBUL-5.5	DE-XI-SCTIUN-3.4-6.5	DE-XI-SCZI-3.4-6.5
Art-No.:	ZFX11000706	ZFX11001128	ZFX09000981	ZFX09000623
<b>ø 6.5 mm</b>	-	S-TI-Base-DE-XI-TBUL-6.5	DE-XI-SCTIUN-3.4-6.5	DE-XI-SCZI-3.4-6.5
Art-No.:	-	ZFX11000707	ZFX09000981	ZFX09000623

TWIN PACK

Components Accessories

	Evolution Matchholder	IntraScan Matchholder H4	IntraScan Matchholder H7	Model analogue
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.0 mm</b>	EVO-DE-XI-MACH-3.0	INTRA-DE-XI-MACH-3.0-H4	INTRA-DE-XI-MACH-3.0-H7	MC-Analog-DE-XI-3.0
Art-No.:	ZFX05002034	ZFX05001974	ZFX05001975	ZFX15002132
<b>ø 3.4 mm</b>	EVO-DE-XI-MACH-3.4	INTRA-DE-XI-MACH-3.4-H4	INTRA-DE-XI-MACH-3.4-H7	MC-Analog-DE-XI-3.4
Art-No.:	ZFX05002035	ZFX05001976	ZFX05001977	ZFX15002133
<b>ø 3.8 mm</b>	EVO-DE-XI-MACH-3.8	INTRA-DE-XI-MACH-3.8-H4	INTRA-DE-XI-MACH-3.8-H7	MC-Analog-DE-XI-3.8
Art-No.:	ZFX05002036	ZFX05001978	ZFX05001979	ZFX15002134
<b>ø 4.5 mm</b>	EVO-DE-XI-MACH-4.5	INTRA-DE-XI-MACH-4.5-H4	INTRA-DE-XI-MACH-4.5-H7	MC-Analog-DE-XI-4.5
Art-No.:	ZFX05002037	ZFX05001980	ZFX05001981	ZFX15002135
<b>ø 5.5 mm</b>	EVO-DE-XI-MACH-5.5	INTRA-DE-XI-MACH-5.5-H4	INTRA-DE-XI-MACH-5.5-H7	MC-Analog-DE-XI-5.5
Art-No.:	ZFX05002038	ZFX05001982	ZFX0500198 3	ZFX15002136
<b>ø 6.5 mm</b>	EVO-DE-XI-MACH-6.5	INTRA-DE-XI-MACH-6.5-H4	INTRA-DE-XI-MACH-6.5-H7	MC-Analog-DE-XI-6.5
Art-No.:	ZFX05002039	ZFX05001984	ZFX05001985	ZFX15002137

Zfx™ Pre-Abutment Blanks including screw




<b>ø 3.0 mm</b>	Art-No.:	ZFX18000260
<b>ø 3.4 mm</b>	Art-No.:	ZFX18000270
<b>ø 3.8 mm</b>	Art-No.:	ZFX18000280
<b>ø 4.5 - 6.5 mm</b>	Art-No.:	ZFX18000290

## ME-ON Type

Compatible with the implant system **AnyOne®** by Megagen.

TWIN PACK

### Components Accessories

	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-			
<b>ø 3.5/7.0 mm</b>	EVO-DU-SU-MACH-3.5-7.0	-	-	-
Art-No.:	ZFX05002202	-	-	-

## ME-AN Type

Compatible with the implant system **AnyRidge®** by Megagen.

### Single Implant Restorations

	<i>Titanium base protected against rotation - Regular Design</i>	<i>Titanium base protected against rotation - Small Design</i>	<i>Screw Standard</i>	<i>Screw for All-Zirconia Abutments</i>
		 € 49.-	 € 19.-	 € 19.-
<b>ø 3.5 / 8.0 mm</b>	-	S-TI-Base-ME-AN-TBLO-3.5-8.0	ME-AN-SCTI-3.5-8.0	ME-AN-SCTI-3.5-8.0
Art-No.:	-	ZFX11002213	ZFX09002258	ZFX09002258
			Recommended Torque 25-35 Ncm	Recommended Torque 25-35 Ncm

TWIN PACK

### Components Accessories

	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.5 / 8.0 mm</b>	EVO-ME-AN-MACH-3.5-8.0	INTRA-ME-AN-MACH-3.5-8.0-H4	INTRA-ME-AN-MACH-3.5-8.0-H7	MC-Analog-ME-AN-3.5-8.0
Art-No.:	ZFX05002269	ZFX05001922	ZFX05001923	ZFX15002326

# ME-EZ Type

Compatible with the implant system **EZ Plus®** by Megagen.

## Single Implant Restorations

	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
		 € 49.-	 € 19.-	 € 19.-
<b>ø 3.3 mm</b>	-	S-TI-Base-ME-EZ-TBLO-3.3	ME-EZ-SCTI-3.3	ME-EZ-SCTI-3.3
Art-No.:	-	ZFX11002188	ZFX09002182	ZFX09002182
<b>ø 4.0 / 5.0 mm</b>	-	S-TI-Base-ME-EZ-TBLO-4.0-5.0	ME-EZ-SCTI-4.0-5.0	ME-EZ-SCTI-4.0-5.0
Art-No.:	-	ZFX11002189	ZFX09002183	ZFX09002183

## Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
		 € 49.-	 € 19.-	 € 19.-
<b>ø 3.3 mm</b>	-	S-TI-Base-ME-EZ-TBUL-3.3	ME-EZ-SCTI-3.3	ME-EZ-SCTI-3.3
Art-No.:	-	ZFX11002191	ZFX09002182	ZFX09002182
<b>ø 4.0 / 5.0 mm</b>	-	S-TI-Base-ME-EZ-TBUL-4.0-5.0	ME-EZ-SCTI-4.0-5.0	ME-EZ-SCTI-4.0-5.0
Art-No.:	-	ZFX11002192	ZFX09002183	ZFX09002183

## Components Accessories

**TWIN PACK**

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-	
<b>ø 3.3 mm</b>	EVO-ME-EZ-MACH-3.3	INTRA-ME-EZ-MACH-3.3-H4	INTRA-ME-EZ-MACH-3.3-H7	MC-Analog-ME-EZ-3.3
Art-No.:	ZFX05002175	ZFX05001232	ZFX05001233	ZFX15002323
<b>ø 4.0 / 5.0 mm</b>	EVO-ME-EZ-MACH-4.0-5.0	INTRA-ME-EZ-MACH-4.0-5.0-H4	INTRA-ME-EZ-MACH-4.0-5.0-H7	MC-Analog-ME-EZ-4.0-5.0
Art-No.:	ZFX05002176	ZFX05001234	ZFX05001235	ZFX15002324

## ME-RE Type

Compatible with the implant system **Rescue Internal®** by Megagen.

### Single Implant Restorations

	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
		€ 49.-	€ 19.-	€ 19.-
<b>ø 4.5 / 6.0 mm</b>	-	S-TI-Base-ME-RE-TBLO-4.5-6.0	ME-RE-SCTI-4.5-6.0	ME-RE-SCTI-4.5-6.0
Art-No.:	-	ZFX11002190	ZFX09002184	ZFX09002184
			25-35 Ncm	25-35 Ncm

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
		€ 49.-	€ 19.-	€ 19.-
<b>ø 4.5 / 6.0 mm</b>	-	S-TI-Base-ME-RE-TBUL-4.5-6.0	ME-RE-SCTI-4.5-6.0	ME-RE-SCTI-4.5-6.0
Art-No.:	-	ZFX11002193	ZFX09002184	ZFX09002184
			25-35 Ncm	25-35 Ncm

### Components Accessories

TWIN PACK

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
2 PC. € 139.-	2 PC. € 149.-	2 PC. € 149.-	2 PC. € 44.-	
<b>ø 4.5 / 6.0 mm</b>	EVO-ME-AN-MACH-3.5-8.0	INTRA-ME-RE-MACH-4.5-6.0-H4	INTRA-ME-RE-MACH-4.5-6.0-H7	MC-Analog-ME-RE-4.5-6.0
Art-No.:	ZFX05002269	ZFX05001100	ZFX05001101	ZFX15002325







## MI-IN Type

Compatible with the implant system **Internal Hex®** by MIS.

**TWIN PACK**






### Components Accessories

	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.3 mm</b>	EVO-MI-IN-MACH-3.3	-	-	-
Art-No.:	ZFX05002204	-	-	-
<b>ø 3.75/4.2 mm</b>	EVO-MI-IN-MACH-3.75-4.2	INTRA-MI-IN-MACH-3.75-4.2-H4	INTRA-MI-IN-MACH-3.75-4.2-H7	MC-Analog-MI-IN-3.75-4.2
Art-No.:	ZFX05002291	ZFX05002293	ZFX05002295	ZFX15002334
<b>ø 5.0/6.0 mm</b>	EVO-MI-IN-MACH-5.0-6.0	INTRA-MI-IN-MACH-5.0-6.0-H4	INTRA-MI-IN-MACH-5.0-6.0-H7	MC-Analog-MI-IN-5.0-6.0
Art-No.:	ZFX05002292	ZFX05002294	ZFX05002296	ZFX15002335






## NO-AC Type

Compatible with the implant system **Active®** by Nobel Biocare.

### Single Implant Restorations





	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
	 € 49.-	 € 49.-	 € 19.-	 € 19.-
<b>ø 3.5 mm</b>	TI-Base-NO-AC-TBLO-3.5	S-TI-Base-NO-AC-TBLO-3.5	NO-AC-SCTI-3.5	NO-AC-SCTI-3.5
Art-No.:	ZFX11000713	ZFX11001130	ZFX09000627	ZFX09000627
<b>ø 4.3 / 5.0 mm</b>	TI-Base-NO-AC-TBLO-4.3-5.0	S-TI-Base-NO-AC-TBLO-4.3-5.0	NO-AC-SCTI-4.3-5.0	NO-AC-SCTI-4.3-5.0
Art-No.:	ZFX11000714	ZFX11001131	ZFX09000628	ZFX09000628
			<b>Recommended Torque</b> 35 Ncm	<b>Recommended Torque</b> 35 Ncm

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
	 € 49.-	 € 49.-	 € 19.-	 € 19.-
<b>ø 3.5 mm</b>	TI-Base-NO-AC-TBUL-3.5	S-TI-Base-NO-AC-TBUL-3.5	NO-AC-SCTIUN-3.5	NO-AC-SCTIUN-3.5
Art-No.:	ZFX11000715	ZFX11001132	ZFX09000982	ZFX09000982
<b>ø 4.3 / 5.0 mm</b>	TI-Base-NO-AC-TBUL-4.3-5.0	S-TI-Base-NO-AC-TBUL-4.3-5.0	NO-AC-SCTIUN-4.3-5.0	NO-AC-SCTIUN-4.3-5.0
Art-No.:	ZFX11000716	ZFX11001133	ZFX09000983	ZFX09000983
			<b>Recommended Torque</b> 35Ncm	<b>Recommended Torque</b> 35 Ncm

TWIN PACK

Components Accessories

	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.5 mm</b>	EVO-NO-AC-MACH-3.5	INTRA-NO-AC-MACH-3.5-H4	INTRA-NO-AC-MACH-3.5-H7	MC-Analog-NO-AC-3.5
Art-No.:	ZFX05002040	ZFX05001986	ZFX05001987	ZFX15002138
<b>ø 4.3 / 5.0 mm</b>	EVO-NO-AC-MACH-4.3-5.0	INTRA-NO-AC-MACH-4.3-5.0-H4	INTRA-NO-AC-MACH-4.3-5.0-H7	MC-Analog-NO-AC-4.3-5.0
Art-No.:	ZFX05002041	ZFX05001988	ZFX05001989	ZFX15002139

Zfx™ Pre-Abutment Blanks including screw



<b>ø 3.5 mm</b>	Art-No.:	ZFX18000210
<b>ø 4.3 - 5.0 mm</b>	Art-No.:	ZFX18000220

## NO-BR Type

Compatible with the implant system **Brånemark®** by Nobel Biocare.

### Single Implant Restorations





	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
	 € 49.-	 € 49.-	 € 19.-	 € 19.-
<b>∅ 3.5 mm</b>	TI-Base-NO-BR-TBLO-3.5	S-TI-Base-NO-BR-TBLO-3.5	NO-BR-SCTI-3.5	NO-AC-SCTI-3.5
Art-No.:	ZFX11000717	ZFX11001134	ZFX09000631	ZFX09000627
<b>∅ 4.1 mm</b>	TI-Base-NO-BR-TBLO-4.1	S-TI-Base-NO-BR-TBLO-4.1	NO-BR-SCTI-4.1	NO-AC-SCTI-4.3-5.0
Art-No.:	ZFX11000718	ZFX11001135	ZFX09000632	ZFX09000628
<b>∅ 5.1 mm</b>	TI-Base-NO-BR-TBLO-5.1	S-TI-Base-NO-BR-TBLO-5.1	NO-BR-SCTI-5.1	NO-BR-SCZI-5.1
Art-No.:	ZFX11000719	ZFX11001136	ZFX09000633	ZFX09000636

### Implant Bar / Bridge Restorations


	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
		 € 49.-	 € 19.-	 € 19.-
<b>∅ 3.5 mm</b>	-	S-TI-Base-NO-BR-TBUL-3.5	NO-BR-SCTIUN-3.5	NO-AC-SCTI-3.5
Art-No.:	-	ZFX11000720	ZFX09001166	ZFX09000627
<b>∅ 4.1 mm</b>	-	S-TI-Base-NO-BR-TBUL-4.1	NO-BR-SCTI-4.1	NO-AC-SCTI-4.3-5.0
Art-No.:	-	ZFX11000721	ZFX09000632	ZFX09000628
<b>∅ 5.1 mm</b>	-	S-TI-Base-NO-BR-TBUL-5.1	NO-BR-SCTI-5.1	NO-BR-SCZI-5.1
Art-No.:	-	ZFX11000722	ZFX09000633	ZFX09000636

TWIN PACK

Components Accessories

	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139,-	 2 PC. € 149,-	 2 PC. € 149,-	 2 PC. € 44,-
<b>ø 3.5 mm</b>	EVO-NO-BR-MACH-3.5	INTRA-NO-BR-MACH-3.5-H4	INTRA-NO-BR-MACH-3.5-H7	MC-Analog-NO-BR-3.5
Art-No.:	ZFX05002042	ZFX05001990	ZFX05001991	ZFX15002140
<b>ø 4.1 mm</b>	EVO-NO-BR-MACH-4.1	INTRA-NO-BR-MACH-4.1-H4	INTRA-NO-BR-MACH-4.1-H7	MC-Analog-NO-BR-4.1
Art-No.:	ZFX05002043	ZFX05001992	ZFX05001993	ZFX15002141
<b>ø 5.1 mm</b>	EVO-NO-BR-MACH-5.1	INTRA-NO-BR-MACH-5.1-H4	INTRA-NO-BR-MACH-5.1-H7	MC-Analog-NO-BR-5.1
Art-No.:	ZFX05002044	ZFX05001994	ZFX05001995	ZFX15002142

Zfx™ Pre-Abutment Blanks including screw

 € 44.90	<b>ø 3.5 mm</b>	Art-Nr.:	ZFX18000230
	<b>ø 4.1 mm</b>	Art-Nr.:	ZFX18000240
	<b>ø 5.1 mm</b>	Art-Nr.:	ZFX18000250

## NO-MU Type

Compatible with the implant system **Multi-unit®** by Nobel Biocare.

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
	 € 49.-		 € 19.-	 € 19.-
<b>ø 4.8 mm</b>	TI-Base-NO-MU-TBUL-4.8	-	NO-MU-SCTI-4.8	NO-MU-SCTI-4.8
Art-No.:	ZFX11000723	-	ZFX09000637	ZFX09000637
<b>ø 6.0 mm</b>	TI-Base-NO-MU-TBUL-6.0	-	NO-MU-SCTI-6.0	NO-MU-SCTI-6.0
Art-No.:	ZFX11000724	-	ZFX09000638	ZFX09000638

### Components Accessories TWIN PACK

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
 2 PC. € 139.-	 2 PC. € 149.-		 2 PC. € 44.-	
<b>ø 4.8 mm</b>	EVO-NO-MU-MACH-4.8	INTRA-NO-MU-MACH-4.8-H4	-	MC-Analog-NO-MU-4.8
Art-No.:	ZFX05002045	ZFX05001228	-	ZFX15002300
<b>ø 6.0 mm</b>	EVO-NO-MU-MACH-6.0	INTRA-NO-MU-MACH-6.0-H4	-	MC-Analog-NO-MU-6.0
Art-No.:	ZFX05002046	ZFX05001230	-	ZFX15002301

# Uncompromising Perfection

Innovative CAD/CAM technologies  
for dental practices and dental laboratories

★ ★ ★  
THE DENTAL  
SCAN CHAMPION  
2017



The **Zfx™ Evolution plus** is a high-end dental scanner of the latest generation and due to the absolute measuring accuracy of less than 9 µm in the solid (140 x 80 mm) according to VDI test procedures, it's one of the most accurate dental scanner on the market. Equipped with the most modern technologies and the knowledge after decades of digital pioneer work, it enables best fitting for directly screwed works on

implant level! In the field of 5-axis simultaneous machining the **Zfx™ Inhouse5x**, consisting of a steel body weighing 220kg, guarantees maximum stability and highest accuracy in a tiny space. The milling unit is suitable for dry and wet machining and thus enables also the processing of hard materials such as titanium, cobalt chrome and glass ceramics.

More Information on [www.zfx-dental.com](http://www.zfx-dental.com)

The Zfx logo consists of the letters 'Zfx' in a bold, blue, sans-serif font. The 'Z' and 'f' are connected, and the 'x' is separate. The 'f' has a long, thin vertical stroke that extends downwards.

## NO-RE Type

Compatible with the implant system **Replace®** by Nobel Biocare.

### Single Implant Restorations

	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
	 € 49.-	 € 49.-	 € 19.-	 € 19.-
<b>ø 3.5 mm</b>	TI-Base-NO-RE-TBLO-3.5	S-TI-Base-NO-RE-TBLO-3.5	NO-RE-SCTI-3.5	NO-RE-SCZI-3.5
Art-No.:	ZFX11000725	ZFX11001140	ZFX09000641	ZFX09000643
<b>ø 4.3 mm</b>	TI-Base-NO-RE-TBLO-4.3	S-TI-Base-NO-RE-TBLO-4.3	NO-RE-SCTI-4.3-6.0	NO-RE-SCZI-4.3-6.0
Art-No.:	ZFX11000726	ZFX11001141	ZFX09000642	ZFX09000644
<b>ø 5.0 mm</b>	TI-Base-NO-RE-TBLO-5.0	S-TI-Base-NO-RE-TBLO-5.0	NO-RE-SCTI-4.3-6.0	NO-RE-SCZI-4.3-6.0
Art-No.:	ZFX11000727	ZFX11001142	ZFX09000642	ZFX09000644
<b>ø 6.0 mm</b>	TI-Base-NO-RE-TBLO-6.0	S-TI-Base-NO-RE-TBLO-6.0	NO-RE-SCTI-4.3-6.0	NO-RE-SCZI-4.3-6.0
Art-No.:	ZFX11000728	ZFX11001143	ZFX09000642	ZFX09000644





### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
	 € 49.-	 € 49.-	 € 19.-	 € 19.-
<b>ø 3.5 mm</b>	TI-Base-NO-RE-TBUL-3.5	S-TI-Base-NO-RE-TBUL-3.5	NO-RE-SCTI-3.5	NO-RE-SCZI-3.5
Art-No.:	ZFX11000729	ZFX11001144	ZFX09000641	ZFX09000643
<b>ø 4.3 mm</b>	TI-Base-NO-RE-TBUL-4.3	S-TI-Base-NO-RE-TBUL-4.3	NO-RE-SCTI-4.3-6.0	NO-RE-SCZI-4.3-6.0
Art-No.:	ZFX11000730	ZFX11001145	ZFX09000642	ZFX09000644
<b>ø 5.0 mm</b>	TI-Base-NO-RE-TBUL-5.0	S-TI-Base-NO-RE-TBUL-5.0	NO-RE-SCTI-4.3-6.0	NO-RE-SCZI-4.3-6.0
Art-No.:	ZFX11000731	ZFX11001146	ZFX09000642	ZFX09000644
<b>ø 6.0 mm</b>	TI-Base-NO-RE-TBUL-6.0	S-TI-Base-NO-RE-TBUL-6.0	NO-RE-SCTI-4.3-6.0	NO-RE-SCZI-4.3-6.0
Art-No.:	ZFX11000732	ZFX11001147	ZFX09000642	ZFX09000644




TWIN PACK

Components Accessories

	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.5 mm</b>	EVO-NO-RE-MACH-3.5	INTRA-NO-RE-MACH-3.5-H4	INTRA-NO-RE-MACH-3.5-H7	MC-Analog-NO-RE-3.5
Art-No.:	ZFX05002047	ZFX05001996	ZFX05001997	ZFX15002143
<b>ø 4.3 mm</b>	EVO-NO-RE-MACH-4.3	INTRA-NO-RE-MACH-4.3-H4	INTRA-NO-RE-MACH-4.3-H7	MC-Analog-NO-RE-4.3
Art-No.:	ZFX05002048	ZFX05001998	ZFX05001999	ZFX15002144
<b>ø 5.0 mm</b>	EVO-NO-RE-MACH-5.0	INTRA-NO-RE-MACH-5.0-H4	INTRA-NO-RE-MACH-5.0-H7	MC-Analog-NO-RE-5.0
Art-No.:	ZFX05002049	ZFX05002000	ZFX05002001	ZFX15002145
<b>ø 6.0 mm</b>	EVO-NO-RE-MACH-6.0	INTRA-NO-RE-MACH-6.0-H4	INTRA-NO-RE-MACH-6.0-H7	MC-Analog-NO-RE-6.0
Art-No.:	ZFX05002050	ZFX05002002	ZFX05002003	ZFX15002146

Zfx™ Pre-Abutment Blanks including screw

 € 44.90	<b>ø 3.5 mm</b>	Art-No.:	ZFX18000170
	<b>ø 4.3 mm</b>	Art-No.:	ZFX18000180
	<b>ø 5.0 mm</b>	Art-No.:	ZFX18000190
	<b>ø 6.0 mm</b>	Art-No.:	ZFX18000200

## OS-GS Type

Compatible with the implant system **GS & TS System®** by Osstem.

### Single Implant Restorations

	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
<b>ø 3.5 mm</b>	-	S-TI-Base-OS-GS-TBLO-3.5	OS-GS-SCTI-3.5	OS-GS-SCTI-3.5
Art-No.:	-	ZFX11002112	ZFX09002110	ZFX09002110
<b>ø 4.0 / 7.0 mm</b>	-	S-TI-Base-OS-GS-TBLO-4.0-7.0	OS-GS-SCTI-4.0-7.0	OS-GS-SCTI-4.0-7.0
Art-No.:	-	ZFX11002113	ZFX09002111	ZFX09002111

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
<b>ø 3.5 mm</b>	-	S-TI-Base-OS-GS-TBUL-3.5	OS-GS-SCTI-3.5	OS-GS-SCTI-3.5
Art-No.:	-	ZFX11002114	ZFX09002110	ZFX09002110
<b>ø 4.0 / 7.0 mm</b>	-	S-TI-Base-OS-GS-TBUL-4.0-7.0	OS-GS-SCTI-4.0-7.0	OS-GS-SCTI-4.0-7.0
Art-No.:	-	ZFX11002115	ZFX09002111	ZFX09002111

### Components Accessories

**TWIN PACK**

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
<b>ø 3.5 mm</b>	EVO-OS-GS-MACH-3.5	INTRA-OS-GS-MACH-3.5-H4	INTRA-OS-GS-MACH-3.5-H7	MC-Analog-OS-GS-3.5
Art-No.:	ZFX05002173	ZFX05001096	ZFX05001097	ZFX15002308
<b>ø 4.0 / 7.0 mm</b>	EVO-OS-GS-MACH-4.0-7.0	INTRA-OS-GS-MACH-4.0-7.0-H4	INTRA-OS-GS-MACH-4.0-7.0-H7	MC-Analog-OS-GS-4.0-7.0
Art-No.:	ZFX05002174	ZFX05001098	ZFX05001099	ZFX15002309

# PI-AMP Type

Compatible with the implant system **Amplified®** by P-I.

## Single Implant Restorations

	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
			 € 19.-	 € 19.-
			 Recommended Torque	 Recommended Torque
<b>ø 3.5 mm</b>	-	-	PI-AMP-SCTI-3.5-5.1	PI-AMP-SCTI-3.5-5.1
Art-No.:	-	-	101688	101688
<b>ø 4.3 / 5.1 mm</b>	-	-	PI-AMP-SCTI-3.5-5.1	PI-AMP-SCTI-3.5-5.1
Art-No.:	-	-	101688	101688

## Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
			 € 19.-	 € 19.-
			 Recommended Torque	 Recommended Torque
<b>ø 3.5 mm</b>	-	-	PI-AMP-SCTI-3.5-5.1	PI-AMP-SCTI-3.5-5.1
Art-No.:	-	-	101688	101688
<b>ø 4.3 / 5.1 mm</b>	-	-	PI-AMP-SCTI-3.5-5.1	PI-AMP-SCTI-3.5-5.1
Art-No.:	-	-	101688	101688

## Components Accessories

**TWIN PACK**

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.5 mm</b>	EVO-PI-AMP-MACH-3.5	INTRA-PI-AMP-MACH-3.5-H4	INTRA-PI-AMP-MACH-3.5-H7	MC-Analog-PI-AMP-3.5
Art-No.:	ZFX05002273	ZFX05002297	ZFX05002299	ZFX15002336
<b>ø 4.3/5.1 mm</b>	EVO-PI-AMP-MACH-4.3-5.1	INTRA-PI-AMP-MACH-4.3-5.1-H4	INTRA-PI-AMP-MACH-4.3-5.1-H7	MC-Analog-PI-AMP-4.3-5.1
Art-No.:	ZFX05002274	ZFX05002298	ZFX05002300	ZFX15002337

## PI-CO Type

Compatible with the implant system **Conical Abutment® by P-I.**

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
 			 € 19.-	 € 19.-
<b>ø 4.8 mm</b>	-	-	PI-CO-SCTI-4.8	PI-CO-SCTI-4.8
Art-No.:	-	-	101120	101120
			 Recommended Torque 15 Ncm	 Recommended Torque 15 Ncm

### Components Accessories

TWIN PACK

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
	 2 PC. € 139.-	 2 PC. € 149.-		 2 PC. € 44.-
<b>ø 4.8 mm</b>	EVO-PI-CO-MACH-4.8	INTRA-PI-CO-MACH-4.8-H4	-	MC-Analog-PI-CO-4.8
Art-No.:	ZFX05002275	ZFX05002301	-	ZFX15002338

# PI-HE Type

Compatible with the implant system **External Hexagon®** by P-I.

## Single Implant Restorations

	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
			 € 19.-	 € 19.-
<b>ø 3.5 mm</b>	-	-	PI-AMP-SCTI-3.5-5.1   25 Ncm	PI-AMP-SCTI-3.5-5.1   25 Ncm
Art-No.:	-	-	101688	101688
<b>ø 4.1 mm</b>	-	-	PI-HE-SCTI-4.1-5.1   35 Ncm	PI-HE-SCTI-4.1-5.1   35 Ncm
Art-No.:	-	-	101121	101121
<b>ø 5.1 mm</b>	-	-	PI-HE-SCTI-4.1-5.1   35 Ncm	PI-HE-SCTI-4.1-5.1   35 Ncm
Art-No.:	-	-	101121	101121

## Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
			 € 19.-	 € 19.-
<b>ø 3.5 mm</b>	-	-	PI-AMP-SCTI-3.5-5.1   25 Ncm	PI-AMP-SCTI-3.5-5.1   25 Ncm
Art-No.:	-	-	101688	101688
<b>ø 4.1 mm</b>	-	-	PI-HE-SCTI-4.1-5.1   35 Ncm	PI-HE-SCTI-4.1-5.1   35 Ncm
Art-No.:	-	-	101121	101121
<b>ø 5.1 mm</b>	-	-	PI-HE-SCTI-4.1-5.1   35 Ncm	PI-HE-SCTI-4.1-5.1   35 Ncm
Art-No.:	-	-	101121	101121








## Components Accessories TWIN PACK

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-	
<b>ø 3.5 mm</b>	EVO-PI-HE-MACH-3.5	INTRA-PI-HE-MACH-3.5-H4	INTRA-PI-HE-MACH-3.5-H7	MC-Analog-PI-HE-3.5
Art-No.:	ZFX05002276	ZFX05002302	ZFX05002305	ZFX15002339
<b>ø 4.1 mm</b>	EVO-PI-HE-MACH-4.1	INTRA-PI-HE-MACH-4.1-H4	INTRA-PI-HE-MACH-4.1-H7	MC-Analog-PI-HE-4.1
Art-No.:	ZFX05002277	ZFX05002303	ZFX05002306	ZFX15002340
<b>ø 5.1 mm</b>	EVO-PI-HE-MACH-5.1	INTRA-PI-HE-MACH-5.1-H4	INTRA-PI-HE-MACH-5.1-H7	MC-Analog-PI-HE-5.1
Art-No.:	ZFX05002278	ZFX05002304	ZFX05002307	ZFX15002341








## PI-MT Type

Compatible with the implant system **Morse Taper®** by P-I.

### Single Implant Restorations

	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
			 € 22.50.-	 € 22.50.-
<b>ø 3.5 / 5.1 mm</b>	-	-	PI-MT-SCTI-3.5-5.1	PI-MT-SCTI-3.5-5.1
Art-No.:	-	-	171239	171239
			 Recommended Torque 25 Ncm	 Recommended Torque 25 Ncm

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
			 € 22.50.-	 € 22.50.-
<b>ø 3.5 / 5.1 mm</b>	-	-	PI-MT-SCTI-3.5-5.1	PI-MT-SCTI-3.5-5.1
Art-No.:	-	-	171239	171239
			 Recommended Torque 25 Ncm	 Recommended Torque 25 Ncm

**Components Accessories**






**TWIN PACK**

	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 <p>2 PC. € 139.-</p>	 <p>2 PC. € 149.-</p>	 <p>2 PC. € 149.-</p>	 <p>2 PC. € 44.-</p>
<b>ø 3.5 / 5.1 mm</b>	EVO-PI-MT-MACH-3.5-5.1	INTRA-PI-MT-MACH-3.5-5.1-H4	INTRA-PI-MT-MACH-3.5-5.1-H7	MC-Analog-PI-MT-3.5-5.1
Art-No.:	ZFX05002279	ZFX05002308	ZFX05002309	ZFX15002342






## ST-BO Type

Compatible with the implant system **Bone Level®** by Straumann Institut.

### Single Implant Restorations

	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
				
	€ 49.-	€ 49.-	€ 19.-	€ 19.-
<b>∅ 3.3 mm</b>	TI-Base-ST-BO-TBLO-3.3	S-TI-Base-ST-BO-TBLO-3.3	ST-BO-SCTI-3.3	ST-BO-SCZI-3.3
Art-No.:	ZFX11000733	ZFX11001148	ZFX09000645	ZFX09000984
<b>∅ 4.1 / 4.8 mm</b>	TI-Base-ST-BO-TBLO-4.1-4.8	S-TI-Base-ST-BO-TBLO-4.1-4.8	ST-BO-SCTI-4.1-4.8	ST-BO-SCZI-4.1-4.8
Art-No.:	ZFX11000734	ZFX11001149	ZFX09000646	ZFX09000985
			<b>Recommended Torque</b> 35 Ncm	<b>Recommended Torque</b> 35 Ncm





### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
				
	€ 49.-	€ 49.-	€ 19.-	€ 19.-
<b>∅ 3.3 mm</b>	TI-Base-ST-BO-TBUL-3.3	S-TI-Base-ST-BO-TBUL-3.3	ST-BO-SCTI-3.3	ST-BO-SCZI-3.3
Art-No.:	ZFX11000735	ZFX11001150	ZFX09000645	ZFX09000984
<b>∅ 4.1 / 4.8 mm</b>	TI-Base-ST-BO-TBUL-4.1-4.8	S-TI-Base-ST-BO-TBUL-4.1-4.8	ST-BO-SCTI-4.1-4.8	ST-BO-SCZI-4.1-4.8
Art-No.:	ZFX11000736	ZFX11001151	ZFX09000646	ZFX09000985
			<b>Recommended Torque</b> 35 Ncm	<b>Recommended Torque</b> 35 Ncm




Components Accessories

TWIN PACK

	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.3 mm</b>	EVO-ST-BO-MACH-3.3	INTRA-ST-BO-MACH-3.3-H4	INTRA-ST-BO-MACH-3.3-H7	MC-Analog-ST-BO-3.3
Art-No.:	ZFX05002051	ZFX05002004	ZFX05002005	ZFX15002147
<b>ø 4.1 / 4.8 mm</b>	EVO-ST-BO-MACH-4.1-4.8	INTRA-ST-BO-MACH-4.1-4.8-H4	INTRA-ST-BO-MACH-4.1-4.8-H7	MC-Analog-ST-BO-4.1-4.8
Art-No.:	ZFX05002052	ZFX05002006	ZFX05002007	ZFX15002148

Zfx™ Pre-Abutment Blanks including screw

 € 44.90	<b>ø 3.3 mm</b>	Art-No.:	ZFX18000100
	<b>ø 4.1 - 4.8 mm</b>	Art-No.:	ZFX18000110

## ST-MU Type

Compatible with the implant system **Multi-Base®** by Straumann Institut.

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
		€ 49.-	€ 19.-	€ 19.-
<b>ø 3.5 mm</b>	-	S-TI-Base-ST-MU-TBUL-3.5	ST-MU-SCTI-3.5-6.5	15 Ncm
Art-No.:	-	ZFX11002096	ZFX09001011	ZFX09001011
<b>ø 4.5 mm</b>	-	S-TI-Base-ST-MU-TBUL-4.5	ST-MU-SCTI-3.5-6.5	15 Ncm
Art-No.:	-	ZFX11002097	ZFX09001011	ZFX09001011
<b>ø 6.5 mm</b>	-	S-TI-Base-ST-MU-TBUL-6.5	ST-MU-SCTI-3.5-6.5	15 Ncm
Art-No.:	-	ZFX11002098	ZFX09001011	ZFX09001011









### Components Accessories TWIN PACK

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
	2 PC. € 139.-	2 PC. € 149.-		2 PC. € 44.-
<b>ø 3.5 mm</b>	EVO-ST-MU-MACH-3.5	INTRA-ST-MU-MACH-3.5-H4	-	MC-Analog-ST-MU-3.5
Art-No.:	ZFX05002168	ZFX05001934	-	ZFX15002302
<b>ø 4.5 mm</b>	EVO-ST-MU-MACH-4.5	INTRA-ST-MU-MACH-4.5-H4	-	MC-Analog-ST-MU-4.5
Art-No.:	ZFX05002169	ZFX05001935	-	ZFX15002303
<b>ø 6.5 mm</b>	EVO-ST-MU-MACH-6.5	INTRA-ST-MU-MACH-6.5-H4	-	MC-Analog-ST-MU-6.5
Art-No.:	ZFX05002170	ZFX05001936	-	ZFX15002304

# ST-OC Type





Compatible with the implant system **Octa® (External Octagon) by Straumann Institut.**

## Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
 			 € 19.-	 € 19.-
			 Recommended Torque	 Recommended Torque
<b>ø 4.8 mm</b>	-	-	ST-OC-SCTI-4.8-6.5	ST-OC-SCTI-4.8-6.5
Art-No.:	-	-	ZFX09002326	ZFX09002326
<b>ø 6.5 mm</b>	-	-	ST-OC-SCTI-4.8-6.5	ST-OC-SCTI-4.8-6.5
Art-No.:	-	-	ZFX09002326	ZFX09002326

## Components Accessories

**TWIN PACK**

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
 2 PC. € 139.-				
<b>ø 4.8 mm</b>	EVO-ST-OC-MACH-4.8	-	-	-
Art-No.:	ZFX05002280	-	-	-
<b>ø 6.5 mm</b>	EVO-ST-OC-MACH-6.5	-	-	-
Art-No.:	ZFX05002281	-	-	-



## ST-SY Type

Compatible with the implant system **SynOcta®** by Straumann Institut.

### Single Implant Restorations





	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
	 € 49.–	 € 49.–	 € 19.–	 € 19.–
<b>ø 3.5 mm</b>	TI-Base-ST-SY-TBLO-3.5	S-TI-Base-ST-SY-TBLO-3.5	ST-SY-SCTI-3.5	ST-SY-SCZI-3.5
Art-No.:	ZFX11000737	ZFX11001152	ZFX09000648	ZFX09000652
<b>ø 4.8 mm</b>	TI-Base-ST-SY-TBLO-4.8	-	ST-SY-SCTI-4.8-6.5	ST-SY-SCZI-4.8-6.5
Art-No.:	ZFX11000738	-	ZFX09000649	ZFX09000986
<b>ø 6.5 mm</b>	TI-Base-ST-SY-TBLO-6.5	-	ST-SY-SCTI-4.8-6.5	ST-SY-SCZI-4.8-6.5
Art-No.:	ZFX11000739	-	ZFX09000649	ZFX09000986

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
	 € 49.–	 € 49.–	 € 19.–	 € 19.–
<b>ø 3.5 mm</b>	-	S-TI-Base-ST-SY-TBUL-3.5	ST-SY-SCTI-3.5	ST-SY-SCZI-3.5
Art-No.:	-	ZFX11000740	ZFX09000648	ZFX09000652
<b>ø 4.8 mm</b>	TI-Base-ST-SY-TBUL-4.8	-	ST-SY-SCTIUN-4.8-6.5	ST-SY-SCTIUN-4.8-6.5
Art-No.:	ZFX11000741	-	ZFX09000651	ZFX09000651
<b>ø 6.5 mm</b>	TI-Base-ST-SY-TBUL-6.5	-	ST-SY-SCTIUN-4.8-6.5	ST-SY-SCTIUN-4.8-6.5
Art-No.:	ZFX11000742	-	ZFX09000651	ZFX09000651

TWIN PACK

Components Accessories

	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.5 mm</b>	EVO-ST-SY-MACH-3.5	INTRA-ST-SY-MACH-3.5-H4	INTRA-ST-SY-MACH-3.5-H7	MC-Analog-ST-SY-3.5
Art-No.:	ZFX05002053	ZFX05002008	ZFX05002009	ZFX15002149
<b>ø 4.8 mm</b>	EVO-ST-SY-MACH-4.8	INTRA-ST-SY-MACH-4.8-H4	INTRA-ST-SY-MACH-4.8-H7	MC-Analog-ST-SY-4.8
Art-No.:	ZFX05002054	ZFX05002010	ZFX05002011	ZFX15002150
<b>ø 6.5 mm</b>	EVO-ST-SY-MACH-6.5	INTRA-ST-SY-MACH-6.5-H4	INTRA-ST-SY-MACH-6.5-H7	MC-Analog-ST-SY-6.5
Art-No.:	ZFX05002055	ZFX05002012	ZFX05002013	ZFX15002151

Zfx™ Pre-Abutment Blanks including screw







<b>ø 4.8 mm</b>	Art-No.:	ZFX18000120
<b>ø 6.5 mm</b>	Art-No.:	ZFX18000130

## SW-GL Type

Compatible with the implant system **Global® by Sweden & Martina.**

TWIN PACK

### Components Accessories





	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.8 mm</b>	EVO-SW-GL-MACH-3.8	INTRA-SW-GL-MACH-3.8-H4	INTRA-SW-GL-MACH-3.8-H7	MC-Analog-SW-GL-3.8
Art-No.:	ZFX05002261	ZFX05001124	ZFX05001125	ZFX15002315
<b>ø 4.3 mm</b>	EVO-SW-GL-MACH-4.3	INTRA-SW-GL-MACH-4.3-H4	INTRA-SW-GL-MACH-4.3-H7	MC-Analog-SW-GL-4.3
Art-No.:	ZFX05002262	ZFX05001126	ZFX05001127	ZFX15002316
<b>ø 4.8 mm</b>	EVO-SW-GL-MACH-4.8	INTRA-SW-GL-MACH-4.8-H4	INTRA-SW-GL-MACH-4.8-H7	MC-Analog-SW-GL-4.8
Art-No.:	ZFX05002263	ZFX05001128	ZFX05001129	ZFX15002317
<b>ø 5.5 mm</b>	EVO-SW-GL-MACH-5.5	INTRA-SW-GL-MACH-5.5-H4	INTRA-SW-GL-MACH-5.5-H7	MC-Analog-SW-GL-5.5
Art-No.:	ZFX05002264	ZFX05001130	ZFX05001131	ZFX15002318

## SW-OU Type

Compatible with the implant system **Outlink<sup>2</sup>® by Sweden & Martina.**

TWIN PACK

### Components Accessories





	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.3 mm</b>	EVO-SW-OU-MACH-3.3	INTRA-SW-OU-MACH-3.3-H4	INTRA-SW-OU-MACH-3.3-H7	MC-Analog-SW-OU-3.3
Art-No.:	ZFX05002265	ZFX05001142	ZFX05001143	ZFX15002319
<b>ø 4.1 mm SP</b>	EVO-SW-OU-SP-MACH-4.1	INTRA-SW-OU-SP-MACH-4.1-H4	INTRA-SW-OU-SP-MACH-4.1-H7	MC-Analog-SW-OU-SP-4.1
Art-No.:	ZFX05002266	ZFX05001146	ZFX05001147	ZFX15002320
<b>ø 3.75 / 4.1 mm</b>	EVO-SW-OU-MACH-3.75-4.1	INTRA-SW-OU-MACH-3.75-4.1-H4	INTRA-SW-OU-MACH-3.75-4.1-H7	MC-Analog-SW-OU-3.75-4.1
Art-No.:	ZFX05002267	ZFX05001144	ZFX05001145	ZFX15002321
<b>ø 5.0 mm</b>	EVO-SW-OU-MACH-5.0	INTRA-SW-OU-MACH-5.0-H4	INTRA-SW-OU-MACH-5.0-H7	MC-Analog-SW-OU-5.0
Art-No.:	ZFX05002268	ZFX05001148	ZFX05001149	ZFX15002322

## SW-PR Type

Compatible with the implant system **Premium Kohno®** by Sweden & Martina.

**TWIN PACK**

### Components Accessories

	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.3 mm</b>	EVO-SW-PR-MACH-3.3	INTRA-SW-PR-MACH-3.3-H4	INTRA-SW-PR-MACH-3.3-H7	MC-Analog-SW-PR-3.3
Art-No.:	ZFX05002251	ZFX05001114	ZFX05001115	ZFX15002310
<b>ø 3.8 mm</b>	EVO-SW-PR-MACH-3.8	INTRA-SW-PR-MACH-3.8-H4	INTRA-SW-PR-MACH-3.8-H7	MC-Analog-SW-PR-3.8
Art-No.:	ZFX05002252	ZFX05001116	ZFX05001117	ZFX15002311
<b>ø 4.25 mm</b>	EVO-SW-PR-MACH-4.25	INTRA-SW-PR-MACH-4.25-H4	INTRA-SW-PR-MACH-4.25-H7	MC-Analog-SW-PR-4.25
Art-No.:	ZFX05002253	ZFX05001118	ZFX05001119	ZFX15002312
<b>ø 5.0 mm</b>	EVO-SW-PR-MACH-5.0	INTRA-SW-PR-MACH-5.0-H4	INTRA-SW-PR-MACH-5.0-H7	MC-Analog-SW-PR-5.0
Art-No.:	ZFX05002254	ZFX05001120	ZFX05001121	ZFX15002313
<b>ø 6.0 mm</b>	EVO-SW-PR-MACH-6.0	INTRA-SW-PR-MACH-6.0-H4	INTRA-SW-PR-MACH-6.0-H7	MC-Analog-SW-PR-6.0
Art-No.:	ZFX05002255	ZFX05001122	ZFX05001123	ZFX15002314

## TH-SP Type

Compatible with the implant system **SPI® by Thommen Medical.**

### Single Implant Restorations

	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
<b>ø 3.5 mm</b>	-	-	TH-SP-SCTI-3.5   15 Ncm	TH-SP-SCTI-3.5   15 Ncm
Art-No.:	-	1.04.750*	ZFX09002251	ZFX09002251
<b>ø 4.0 mm</b>	-	-	TH-SP-SCTI-4.0-6.0   25 Ncm	TH-SP-SCTI-4.0-6.0   20 Ncm
Art-No.:	-	1.04.751*	ZFX09002250	ZFX09002250
<b>ø 4.5 mm</b>	-	-	TH-SP-SCTI-4.0-6.0   25 Ncm	TH-SP-SCTI-4.0-6.0   20 Ncm
Art-No.:	-	1.04.752*	ZFX09002250	ZFX09002250
<b>ø 5.0 mm</b>	-	-	TH-SP-SCTI-4.0-6.0   25 Ncm	TH-SP-SCTI-4.0-6.0   20 Ncm
Art-No.:	-	1.04.753*	ZFX09002250	ZFX09002250
<b>ø 6.0 mm</b>	-	-	TH-SP-SCTI-4.0-6.0   25 Ncm	TH-SP-SCTI-4.0-6.0   20 Ncm
Art-No.:	-	1.04.754*	ZFX09002250	ZFX09002250

### Implant Bar / Bridge Restorations





	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b>	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
<b>ø 3.5 mm</b>	-	-	TH-SP-SCTI-3.5   15 Ncm	TH-SP-SCTI-3.5   15 Ncm
Art-No.:	-	1.04.760*	ZFX09002251	ZFX09002251
<b>ø 4.0 mm</b>	-	-	TH-SP-SCTI-4.0-6.0   25 Ncm	TH-SP-SCTI-4.0-6.0   20 Ncm
Art-No.:	-	1.04.761*	ZFX09002250	ZFX09002250
<b>ø 4.5 mm</b>	-	-	TH-SP-SCTI-4.0-6.0   25 Ncm	TH-SP-SCTI-4.0-6.0   20 Ncm
Art-No.:	-	1.04.762*	ZFX09002250	ZFX09002250
<b>ø 5.0 mm</b>	-	-	TH-SP-SCTI-4.0-6.0   25 Ncm	TH-SP-SCTI-4.0-6.0   20 Ncm
Art-No.:	-	1.04.763*	ZFX09002250	ZFX09002250
<b>ø 6.0 mm</b>	-	-	TH-SP-SCTI-4.0-6.0   25 Ncm	TH-SP-SCTI-4.0-6.0   20 Ncm
Art-No.:	-	1.04.764*	ZFX09002250	ZFX09002250

\*Not available at Zfx. This code is the original Thommen Medical's article code.



## Components Accessories

TWIN PACK


	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.5 mm</b>	EVO-TH-SP-MACH-3.5	INTRA-TH-SP-MACH-3.5-H4	INTRA-TH-SP-MACH-3.5-H7	MC-Analog-TH-SP-3.5
Art-No.:	ZFX05002246	ZFX05001104	ZFX05001105	ZFX15002284
<b>ø 4.0 mm</b>	EVO-TH-SP-MACH-4.0	INTRA-TH-SP-MACH-4.0-H4	INTRA-TH-SP-MACH-4.0-H7	MC-Analog-TH-SP-4.0
Art-No.:	ZFX05002247	ZFX05001106	ZFX05001107	ZFX15002285
<b>ø 4.5 mm</b>	EVO-TH-SP-MACH-4.5	INTRA-TH-SP-MACH-4.5-H4	INTRA-TH-SP-MACH-4.5-H7	MC-Analog-TH-SP-4.5
Art-No.:	ZFX05002248	ZFX05001108	ZFX05001109	ZFX15002286
<b>ø 5.0 mm</b>	EVO-TH-SP-MACH-5.0	INTRA-TH-SP-MACH-5.0-H4	INTRA-TH-SP-MACH-5.0-H7	MC-Analog-TH-SP-5.0
Art-No.:	ZFX05002249	ZFX05001110	ZFX05001111	ZFX15002287
<b>ø 6.0 mm</b>	EVO-TH-SP-MACH-6.0	INTRA-TH-SP-MACH-6.0-H4	INTRA-TH-SP-MACH-6.0-H7	MC-Analog-TH-SP-6.0
Art-No.:	ZFX05002250	ZFX05001112	ZFX05001113	ZFX15002288

## ZI-EZ Type


Compatible with the implant system **Eztetic Implant®** by Zimmer.

## Components Accessories

TWIN PACK

	<i>Evolution Matchholder</i>	<i>IntraScan Matchholder H4</i>	<i>IntraScan Matchholder H7</i>	<i>Model analogue</i>
	 2 PC. € 139.-			
<b>ø 2.9 mm</b>	EVO-ZI-EZ-MACH-2.9	-	-	-
Art-No.:	ZFX05002270	-	-	-

## Zfx™ Pre-Abutment Blanks including screw

	<b>ø 2.9 mm</b>	Art-No.:	ZFX18000420
€ 44.90			

## ZI-SC Type

Compatible with the implant system **Tapered Screw Vent®** by Zimmer.

### Single Implant Restorations





	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
	 € 49.-	 € 49.-	 € 19.-	 € 19.-
<b>ø 3.5 mm</b>	TI-Base-ZI-SC-TBLO-3.5	S-TI-Base-ZI-SC-TBLO-3.5	ZI-SC-SCTI-3.5-5.7	Screw AH20S
Art-No.:	ZFX11000745	ZFX11001154	ZFX09000655	AH20S
<b>ø 4.5 mm</b>	TI-Base-ZI-SC-TBLO-4.5	S-TI-Base-ZI-SC-TBLO-4.5	ZI-SC-SCTI-3.5-5.7	Screw AH20S
Art-No.:	ZFX11000746	ZFX11001155	ZFX09000655	AH20S
<b>ø 5.7 mm</b>	TI-Base-ZI-SC-TBLO-5.7	S-TI-Base-ZI-SC-TBLO-5.7	ZI-SC-SCTI-3.5-5.7	Screw AH20S
Art-No.:	ZFX11000747	ZFX11001156	ZFX09000655	AH20S

### Implant Bar / Bridge Restorations


	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
	 € 49.-	 € 49.-	 € 19.-	 € 19.-
<b>ø 3.5 mm</b>	TI-Base-ZI-SC-TBUL-3.5	S-TI-Base-ZI-SC-TBUL-3.5	Screw AH20S	Screw AH20S
Art-No.:	ZFX11000748	ZFX11001157	AH20S	AH20S
<b>ø 4.5 mm</b>	TI-Base-ZI-SC-TBUL-4.5	S-TI-Base-ZI-SC-TBUL-4.5	Screw AH20S	Screw AH20S
Art-No.:	ZFX11000749	ZFX11001158	AH20S	AH20S
<b>ø 5.7 mm</b>	TI-Base-ZI-SC-TBUL-5.7	S-TI-Base-ZI-SC-TBUL-5.7	Screw AH20S	Screw AH20S
Art-No.:	ZFX11000750	ZFX11001159	AH20S	AH20S

Components Accessories

TWIN PACK

	Evolution Matchholder	IntraScan Matchholder H4	IntraScan Matchholder H7	Model analogue
	 2 PC. € 139.-	 2 PC. € 149.-	 2 PC. € 149.-	 2 PC. € 44.-
<b>ø 3.5 mm</b>	EVO-ZI-SC-MACH-3.5 Art-No.: ZFX05002056	INTRA-ZI-SC-MACH-3.5-H4 Art-No.: ZFX05001090	INTRA-ZI-SC-MACH-3.5-H7 Art-No.: ZFX05001093	MC-Analog-ZI-SC-3.5 Art-No.: ZFX15002152
<b>ø 4.5 mm</b>	EVO-ZI-SC-MACH-4.5 Art-No.: ZFX05002057	INTRA-ZI-SC-MACH-4.5-H4 Art-No.: ZFX05001091	INTRA-ZI-SC-MACH-4.5-H7 Art-No.: ZFX05001094	MC-Analog-ZI-SC-4.5 Art-No.: ZFX15002153
<b>ø 5.7 mm</b>	EVO-ZI-SC-MACH-5.7 Art-No.: ZFX05002058	INTRA-ZI-SC-MACH-5.7-H4 Art-No.: ZFX05001092	INTRA-ZI-SC-MACH-5.7-H7 Art-No.: ZFX05001095	MC-Analog-ZI-SC-5.7 Art-No.: ZFX15002154

Zfx™ Pre-Abutment Blanks including screw

 € 44.90	<b>ø 2.9 mm</b>	Art-No.:	ZFX18000420
	<b>ø 3.5 mm</b>	Art-No.:	ZFX18000390
	<b>ø 4.5 mm</b>	Art-No.:	ZFX18000400
	<b>ø 5.7 mm</b>	Art-No.:	ZFX18000410

## ZI-SH Type

Compatible with the implant system **Shoulder Abutment®** by Zimmer.

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
		€ 49.-	€ 19.-	€ 19.-
<b>ø 4.5 mm</b>	-	S-TI-Base-ZI-SH-TBUL-4.5	ZI-SH-SCTI-4.5-6.5	ZI-SH-SCTI-4.5-6.5
Art-No.:	-	ZFX11002084	ZFX09001006	ZFX09001006
<b>ø 6.5 mm</b>	-	S-TI-Base-ZI-SH-TBUL-6.5	ZI-SH-SCTI-4.5-6.5	ZI-SH-SCTI-4.5-6.5
Art-No.:	-	ZFX11002085	ZFX09001006	ZFX09001006

### Components Accessories






TWIN PACK

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
	2 PC. € 139.-	2 PC. € 149.-		2 PC. € 44.-
<b>ø 4.5 mm</b>	EVO-ZI-SH-MACH-4.5	INTRA-ZI-SH-MACH-4.5-H4	-	MC-Analog-ZI-SH-4.5
Art-No.:	ZFX05002171	ZFX05001937	-	ZFX15002306
<b>ø 6.5 mm</b>	EVO-ZI-SH-MACH-6.5	INTRA-ZI-SH-MACH-6.5-H4	-	MC-Analog-ZI-SH-6.5
Art-No.:	ZFX05002172	ZFX05001938	-	ZFX15002307






# ZI-SW Type

Compatible with the implant system **SwissPlus®** by Zimmer.

## Single Implant Restorations





	<b>Titanium base protected against rotation</b> - Regular Design	<b>Titanium base protected against rotation</b> - Small Design	<b>Screw</b> Standard	<b>Screw</b> for All-Zirconia Abutments
				
		€ 49.-	€ 19.-	€ 19.-
			$\vec{T}$ Recommended Torque	$\vec{T}$ Recommended Torque
<b>ø 3.8 mm</b>	-	S-TI-Base-ZI-SW-TBLO-3.8	ZI-SW-SCTI-3.8	ZI-SW-SCZI-3.8
Art-No.:	-	ZFX11001160	ZFX09000987	ZFX09001004
<b>ø 4.8 mm</b>	-	S-TI-Base-ZI-SW-TBLO-4.8	ZI-SW-SCTI-4.8	ZI-SW-SCZI-4.8
Art-No.:	-	ZFX11001161	ZFX09000988	ZFX09001005

## Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
				
		€ 49.-	€ 19.-	€ 19.-
			$\vec{T}$ Recommended Torque	$\vec{T}$ Recommended Torque
<b>ø 3.8 mm</b>	-	S-TI-Base-ZI-SW-TBUL-3.8	ZI-SW-SCZI-3.8	ZI-SW-SCZI-3.8
Art-No.:	-	ZFX11001162	ZFX09001004	ZFX09001004
<b>ø 4.8 mm</b>	-	S-TI-Base-ZI-SW-TBUL-4.8	ZI-SW-SCZI-4.8	ZI-SW-SCZI-4.8
Art-No.:	-	ZFX11001163	ZFX09001005	ZFX09001005

## Components Accessories

**TWIN PACK**

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
				
2 PC. € 139.-	2 PC. € 149.-	2 PC. € 149.-	2 PC. € 44.-	
<b>ø 3.8 mm</b>	EVO-ZI-SW-MACH-3.8	INTRA-ZI-SW-MACH-3.8-H4	INTRA-ZI-SW-MACH-3.8-H7	MC-Analog-ZI-SW-3.8
Art-No.:	ZFX05002059	ZFX05002014	ZFX05002015	ZFX15002155
<b>ø 4.8 mm</b>	EVO-ZI-SW-MACH-4.8	INTRA-ZI-SW-MACH-4.8-H4	INTRA-ZI-SW-MACH-4.8-H7	MC-Analog-ZI-SW-4.8
Art-No.:	ZFX05002060	ZFX05002016	ZFX05002017	ZFX15002156

## ZI-TA Type

Compatible with the implant system **Tapered Abutment® by Zimmer.**

### Implant Bar / Bridge Restorations

	<b>Titanium base not protected against rotation</b> - Regular Design	<b>Titanium base not protected against rotation</b> - Small Design	<b>Bridge screw</b> Standard	<b>Bridge screw</b> for All-Zirconia Implant Bridges
		€ 49.-	€ 19.-	€ 19.-
<b>ø 4.5 mm</b>	-	S-TI-Base-ZI-TA-TBUL-4.5	ZI-TA-SCTI-4.5	ZI-TA-SCTI-4.5
Art-No.:	-	ZFX11002078	ZFX09000657	ZFX09000657
			<b>Recommended Torque</b> 20 Ncm	<b>Recommended Torque</b> 20 Ncm

### Components Accessories TWIN PACK

	<b>Evolution Matchholder</b>	<b>IntraScan Matchholder H4</b>	<b>IntraScan Matchholder H7</b>	<b>Model analogue</b>
	2 PC. € 139.-	2 PC. € 149.-		2 PC. € 44.-
<b>ø 4.5 mm</b>	EVO-ZI-TA-MACH-4.5	INTRA-ZI-TA-MACH-4.5-H4	-	MC-Analog-ZI-TA-4.5
Art-No.:	ZFX05002157	ZFX05001924	-	ZFX15002305

## Torque wrench for laboratories

**Torque wrench**  
for Zfx™ Evolution Matchholders



Art-No.: ZFX02002067

## Hex drivers for dentists

**Hex driver, Long Standard**  
for Zfx™ IntraScan Matchholders



Art-No.: HXL1.25

Torque wrench adapter without Gemlock

**Hex driver, Long with Gemlock Retention**  
for Zfx™ IntraScan Matchholders



Art-No.: HXLGR1.25

with Gemlock Retention

# Instructions for use & Warranties

## Bonding technique with titanium base (link)

- ✗ The titanium base acts as an adhesive base for the constructed individual abutment
- ✗ The position of the titanium base is scanned with the corresponding matchholder
- ✗ The scanning process is carried out with the optical scanner Zfx™ Evolution plus
- ✗ The recorded geometry is used for crown construction within the Zfx™ CAD software

## Indications for the titanium base

- ✗ For individual abutments for bonding (e.g. zirconium abutment)
- ✗ With the titanium base "non-engaging", implant bridge constructions (implant bridges) are also feasible

## Indications for matchholder (scan body)

- ✗ To determine the exact position of the implant on the model or the patient
- ✗ To use with the Torque Wrench for Matchholders  
Article No.: ZFX02002067

## Contra-Indications

- ✗ Abutments with an incline correction of over 20° with respect to the implant axis
- ✗ Customized abutments on titanium base with a different diameter to the implant base are not allowed, as this can cause irritation of the peri-implant soft tissue
- ✗ The torque of each screw must match with the size indications stated by the respective implant manufacturer (e.g. ZI-SC-Type = Zimmer Biomet)
- ✗ Implant bridge constructions (implant bridges) with rotation-locked geometry

# Warranties

## General terms and conditions for parts of the implant components

In regard with the supplied parts (titanium base, abutment screw and scan body) the mechanical stability is guaranteed (guarantee secured by legal law). This guarantee is only coming into effect in case that all above-mentioned technical/clinical indications have been absolutely strictly observed. These parts are subject to continuous further development. Therefore we reserve the right to carry out any product modifications without prior notice.

## Accessories

### Laboratory screw and scan base

The enclosed bolt can be used for definitive fixation inside the mouth. For laboratory work it's advisable to use the compatible laboratory screw. For the fixation of the screw please always use the screwdriver of the implant manufacturer.

## Materials

- ✗ Titanium base: titanium alloy (Ti6Al4V)
- ✗ Implant Screw: titanium alloy (Ti6Al4V)
- ✗ Scan body: polyether-ether ketones (PEEK)

## Indications for the Matchholder IntraScan

- ✗ To determine the exact implant position in the oral cavity
- ✗ To use with the Zimmer Hex driver  
Article No.: HXL1.25 or HXLGR1.25

## Indications for the Model analogues

- ✗ The model analogues have to be used for the models, which have been digitally created with the Zfx™ Modelcreator



## Using the matchholders (scan body)

- ✗ To record the implant position during the scanning process, use a scan body, also named as matchholder
- ✗ This must be selected according to the abutment series being applied, compatible with the original implant system and implant diameter
- ✗ The correct position of the matchholders corresponds to the exact position of the implant
- ✗ As soon as positioning is correct, there's no gap existing between matchholder and implant analogue and thus no longer rotation is feasible
- ✗ The matchholder is fixated by the abutment screw included within the delivery

## Processing of titanium base

- ✗ The inner side of the titanium base may not be treated mechanically or with a sandblaster
- ✗ The diameter and length of the titanium base must not be reduced
- ✗ The surface which comes into contact with the zirconium oxide abutment is blasted with aluminium oxide of 2 bar (100 µ) and then cleaned with a steam jet
- ✗ To protect the inside of the titanium base during the sand blasting, it is advisable to fix it on a laboratory analogue
- ✗ The surface conditioning of the titanium base is carried out in accordance with the manufacturer instructions of the respective used adhesive system (e.g. Rocotec®)

## Gluing the components

- ✗ To glue the titanium base with the individual abutments, e.g. "PANAVIA™ F 2.0" (made by Kuraray Europe), RelayXUnicem® (3M-Espe) or other similar systems a suitable compound system is used
- ✗ It is recommended to use the related metal primer
- ✗ The titanium base is fixed on the laboratory analogue by using an abutment screw. The screw channel is sealed with wax **1**
- ✗ Prepare adhesive material and put a thin coat on the titanium base
- ✗ The individual abutments is pushed onto titanium base until the abutment and titanium base are in line with the bearing surfaces **2**
- ✗ Ensure that the abutment is pushed into its correct position
- ✗ Remove the large surpluses of the compound system
- ✗ If required, a plastic tool for fitting the abutment can be created

## Preparation and polishing

- ✗ Once the adhesive has hardened, remove excess with silicone polishers. It's of essential importance that only a thin layer of adhesive remains between the parts. The transition is highly polished by a silicone polishers for ceramic. **3**
- ✗ Clean the screw channel. **4**



# Order form



**Delivery:**  Standard  Express (+ € 29.00 Service fee)

**Yes, I want to benefit from the € 100.- discount.** Minimum order value: € 1,000.-\*

**Firma:**

---

No.	Zfx Art-No.	Description	Quantity	Price in Euro
1.				
2.				
3.				
4.				
5.				
6.				
7.				
8.				
9.				
10.				
11.				
12.				
13.				
14.				
15.				

**Fax +49 8131 33 244 10 or via e-mail: [order@zfx-dental.com](mailto:order@zfx-dental.com)**

Date:

---

Signature:

---

Company stamp

\*There is only one coupon redeemable per order.

# Zfx™ Process Accuracy

Highest fitting without sending models and black box!



The two-piece Zfx™ Scan bodies for multiple usage without any deficit of accuracy are produced with an innovative code system, that enables the automatic recognition of the implant type and prevents application errors. Thus, in combination with the own-developed torque wrench for the Zfx™ Evolution plus, an accuracy is ensured of  $< 5 \mu\text{m}$

The Zfx™ Evolution plus scanner enables a measuring accuracy\* in the solid of  $< 9 \mu\text{m}$

PROCESS  
ACCURACY  
 $< 20 \mu\text{m}$



The industrial manufacturing with 5-axis simultaneous machining achieves an accuracy of  $< 5 \mu\text{m}$

The Zfx logo, with 'Zfx' in a stylized blue font and 'fx' in a black font.

\* Measurement according to the VDI test procedure



[www.zfx-dental.com](http://www.zfx-dental.com)

**Zfx GmbH**

Kopernikusstraße 27  
85221 Dachau, Germany  
Tel. +49 81 31 33 244 - 0  
Fax +49 81 31 33 244 - 10  
[office@zfx-dental.com](mailto:office@zfx-dental.com)